

# SenSmart™ Model 8204CA

*Dual Emitter Regional Oximetry  
(rSO<sub>2</sub>) Sensor*



## Greater Accuracy and Reliability for Enhanced Patient Outcomes


An advancement in near-infrared spectroscopy (NIRS) regional oximetry, the SenSmart Model 8204CA dual emitter regional oximetry sensor for use with the SenSmart Universal Oximetry (rSO<sub>2</sub>/SpO<sub>2</sub>) System, offers clinicians absolute accuracy to optimize patient care. SenSmart Technology offers accurate rSO<sub>2</sub> values on all patients...no more guessing when the system doesn't provide a reading, and no more outliers that misrepresent patients' true physiologic condition.<sup>1,2</sup>

# SenSmart™ rSO<sub>2</sub> Technology: Four-Wavelength Optics


## Greater Accuracy and Reliability through Better Science

SenSmart rSO<sub>2</sub> Technology utilizes four wavelengths of light to measure the balance of oxy- and deoxy-hemoglobin, while compensating for tissue factors that reduce rSO<sub>2</sub> accuracy. The result is highly accurate measurements of patients' true tissue oxygen saturation.

**Conventional Two-Wavelength, Single-Emitter NIRS Trending Variability<sup>1</sup>**




**Nonin rSO<sub>2</sub> Technology Four-Wavelength rSO<sub>2</sub> Values Aligned with SavO<sub>2</sub><sup>1</sup>**



*Data rounded up, percentages add up to greater than 100%*



Compared to traditional NIRS trending technology, SenSmart rSO<sub>2</sub> Technology is absolute, providing enhanced accuracy in rSO<sub>2</sub> measurements on all patients, regardless of skin color<sup>1</sup>...no more guessing or over-treatment if patients do not read, or read inaccurately.

**Nonin SenSmart rSO<sub>2</sub> Technology Four-Wavelength Accuracy<sup>1</sup>**



## SenSmart™ rSO<sub>2</sub> Technology: Dual-Emitter Architecture Effectively Removes Surface Tissue Artifacts

Only the patented SenSmart rSO<sub>2</sub> sensors use dual emitters and dual detectors to provide repeatable cerebral cortex measurements that are minimally affected by intervening tissue or surface effects.<sup>3</sup>



A recent independent, peer-reviewed study showed in a statistically significant manner that of the three cerebral oximeters, Nonin's rSO<sub>2</sub> technology had the least interference from extracranial tissue.<sup>3</sup>


Older NIRS technology claimed to be absolute, but clinician confidence suffered when the technology could not read some patients. In a recent series of 50 healthy adults, Nonin's rSO<sub>2</sub> technology captured data on all subjects, regardless of gender, body type or skin color, while 5% of subjects did not read using other technologies.<sup>1</sup>



# SenSmart™ Model X-100 Universal Oximetry (rSO<sub>2</sub>/SpO<sub>2</sub>) System

## Product Highlights

- **Absolute Accuracy**<sup>1</sup> – 4.1 A<sub>rms</sub><sup>\*</sup>
- **Consistent** – Repeatable results for maximum efficiency and time savings<sup>4</sup>
- **Immediate Response** – Near-instant measurements ensure readings from the moment the sensor is placed
- **Real-time Readings** – Sampling rates less than two seconds apart guarantee up-to-date saturation readings



**SenSmart™ X-100**  
Universal Oximetry  
(rSO<sub>2</sub>/SpO<sub>2</sub>) System

*Portable and durable, the SenSmart X-100 Universal Oximetry System has six-channel display and the ability to monitor rSO<sub>2</sub> and SpO<sub>2</sub> at the same time.*

# SenSmart™ Model 8204CA Regional Oximetry (rSO<sub>2</sub>) Sensor Improved rSO<sub>2</sub> Sensor Utility

- 13% smaller footprint than 8004CA
- Quick-connect intermediate cable that can stay connected to signal processor between cases
- Proven, reliable adhesion; even with highly diaphoretic patients
- Nonin dual-emitter absolute accuracy
- Same “over the top” vs side-of-head cable exit angle
- Write-on Sharpie strip on sensor face for time/date of placement (e.g. ICU use)



**INT-100**  
Intermediate Cable

**8204CA**  
Regional Oximetry Sensor

### References:

<sup>1</sup> Data on file at Nonin Medical, Inc. Plymouth, MN.


<sup>2</sup> Fischer G et al; J Thorac Cardiovasc Surg 2011; 141:815.

<sup>3</sup> Davie SN, Grocott HP. Impact of Extracranial Contamination on Regional Cerebral Oxygen Saturation: A Comparison of Three Cerebral Oximetry Technologies. *Anesthesiology*. 2012; 116(4):834–40.

<sup>4</sup> Lobbetael A et al; White Paper 2009; Nonin Medical, Inc.

\*Accuracy for SenSmart™ EQUANOX Advance™ Sensor, Model 8004CA only. Pediatric through adult >40 kg.

Authorized Distributor:



**Nonin Medical, Inc.**  
13700 1st Avenue North  
Plymouth, MN • 55441-5443 • U.S.A.  
**Tel:** +1.763.553.9968 1.800.356.8874  
**Fax:** +1.763.577.5521  
**E-mail:** info@nonin.com

[sensmart.com](http://sensmart.com) [nonin.com](http://nonin.com)

**Nonin Medical B.V.**  
Prins Hendriklaan 26  
1075 BD Amsterdam • Netherlands  
**Tel:** +31 (0)13-79 99 040  
**Fax:** +31 (0)13-79 99 042  
**E-mail:** infointl@nonin.com