



Concentrating on efficiency and convenience with innovative design

- High Ultrafiltration Rate
- No rinse for easy setup
- Low priming volume
- Even blood flow throughout
- A wide range of products to choose from
- High visibility



CAPIOX[®] Hemoconcentrator HC11 / HC05

CAPIOX[®] Hemoconcentrator

Terumo's HC-series has everything needed to concentrate blood with maximum efficiency.

High UFR

Polysulfone membrane material achieves high UFR (Ultrafiltration Rate) for efficient blood concentration in a short time, while restricting removal of albumin, which is vital to the human body.

No-rinse design

No-rinse hollow fibers enable quick set up even during the operation or in case of emergency. In addition, the angle of the blood outlet port is designed to eliminate air bubbles quickly and easily during priming.

Low priming volume

This keeps the HC series compact and light, consequently requiring little storage space and offering easy disposal. The HC05 requires only 35ml priming volume, making it especially beneficial in small volume requirements.

Wide product range

Terumo HC series hemoconcentrators are available in two sizes having different effective membrane areas. Each size is offered with three types of blood ports.

Smooth blood flow

Design innovations in the inlet blood port ensure smooth blood flow and reduce possible stagnation. Also, there is little space between hollow fibers and O-rings around the port. These design innovations ensure smooth, even blood flow throughout, and prevent thrombus formation.



Conventional type

CAPIOX HC series

High visibility

The entire housing, including blood ports, is transparent for easy visual checking.

Sieving Coefficients (In-vitro)

	Vitamin B12	Inulin	Albumin
Solute Molecular Weight (Daltons)	1,350	5,200	66,000
Sieving Coefficients (%)	99.8	97.4	0.6

Inlet conditions: Vitamin B12 and Inulin (Dissolved in purified water),
Flow rate = 200 ml/min, UFR = 40 ml/min,
Albumin: Bovine blood, Hematocrit = 25±1%,
Protein concentration = 60±5 g/L, Temp. = 37°C,
Flow rate = 500 ml/min



Product Description

Reorder Number	Description	Quantity
CX-HC11S16	1.1 m ² Hemoconcentrator with 3/16" slip port	10 pcs./Carton
CX-HC11S	1.1 m ² Hemoconcentrator with 1/4" slip port	10 pcs./Carton
CX-HC11L	1.1 m ² Hemoconcentrator with luer lock	10 pcs./Carton
CX-HC05S16	0.5 m ² Hemoconcentrator with 3/16" slip port	10 pcs./Carton
CX-HC05S	0.5 m ² Hemoconcentrator with 1/4" slip port	10 pcs./Carton
CX-HC05L	0.5 m ² Hemoconcentrator with luer lock	10 pcs./Carton
XX-XH071	Holder for HC11/HC05	-

Specifications

	HC11S			HC05S		
Holder	XX*XH071					
Fiber Material	Polysulfone					
Effective Surface Area	Approx. 1.1 m ²			Approx. 0.5 m ²		
Priming Volume	70 ml			35 ml		
Blood Ports	3/16" Slip	1/4" Slip	Luer Locks	3/16" Slip	1/4" Slip	Luer Locks
Filtrate Port	1/2" (1/4" Adapter)					
Blood Flow Range	Min. 100 - Max. 500 ml/min.*1					
Maximum Transmembrane Pressure	500 mmHg (66.5kPa)					
Method of sterilization	Ethylene Oxide Gas					

*1: At least 50 ml/min is recommended during suspension of ultrafiltration.



Inlet conditions:
Bovine blood, Hematocrit = 25±1%,
Protein concentration = 50±5 g/L,
B.E. = 0±3 mEq/L, Temp. = 37°C

CE 0197



Distributor

TERUMO CORPORATION

44-1, 2-chome, Hatagaya, Shibuya-ku
Tokyo 151-0072, Japan
Phone: 81-3-3374-8111
Fax: 81-3-3374-8196

TERUMO MEDICAL CORPORATION

2100 Cottontail Lane, Somerset
New Jersey 08873, U.S.A.
Phone: 1-732-302-4900
Fax: 1-732-302-3083
Toll free: 1-800-283-7866

TERUMO EUROPE N.V.

Researchpark zone 2 Haasrode
Interleuvenlaan 40, B-3001
Leuven, Belgium
Phone: 32-16-381-211
Fax: 32-16-400-249

®: Registered Trademark
Manufactured by TERUMO CORPORATION in Japan.
© 1998 TERUMO CORPORATION. Printed in Japan. All rights reserved.

