



Advanced Perfusion System 1



System 1

Touch the future...today.

Terumo® Advanced Perfusion System 1


Touch the future...today.

Unprecedented Customization

The sophisticated design of the Terumo System 1 means many things to today's perfusionist: the ability to customize the system's operation with the touch of a finger; an unprecedented degree of flexibility of the system's physical configuration; the ability to upgrade, adapt, and expand the system to meet future needs; a truly intuitive, easy-to-operate user interface; an innovative service application; and, a sleek, compact, modern look.

But most importantly, Terumo System 1's advanced design provides today's perfusionist with the confidence that comes with knowing that Terumo System 1 will be one of the smartest professional investments they'll ever make.


System 1 redefines the notion of customization with quick, simple configurability to easily adapt for new protocols or individual patient needs. You can choose exactly how you're going to manage a case – which devices you'll use, where you'll locate them, how you'll configure the display on the Central Control Monitor. And System 1's software-based design lets you select advanced functions without specialized hardware.



System 1

Customization

- Choose the number and types of devices you need.
- Mount pumps on the system base or on a pole close to the surgical field.
- Rotate pump raceways to optimize tubing lengths.
- Control pumps using the front panel knob or from the Central Control Monitor.
- Select from expanded choices for pump responses to alert and alarm conditions.
- View pump information either on the front panel display or on the Central Control Monitor.
- Create up to 12 procedure-specific display screens.
- Add devices and sensors as your needs change.



System Modules Plug In Anywhere



System 1's safety and monitoring systems are housed in individual compact modules that can be plugged into any of 18 slots conveniently housed in the system's base. This modular system not only enhances the system's configurability, it also provides an uncluttered and less-distracting working environment.



Select as many as 18 modules that monitor or operate up to:

- 8 pumps (up to 2 centrifugal pumps)
- 1 electronic O₂ blender/analyzer
- 1 occluder
- 4 air bubble detectors
- 1 level alert and 1 level alarm sensor
- 8 pressure sensors
- 8 temperature sensors
- 4 flow sensors
- 1 serial interface that connects to either the CDI™ System 100 or CDI™ System 500
- 1 serial interface that connects directly to the Sarns™ CCPro Data™ Perfusion Software

Compact Base

- Optional electronic O₂ blender includes a built-in analyzer.
- Integrated battery backup can power the system without interruption for one hour.
- Side and back covers protect cables and modules from spills.
- Large, smooth surfaces are designed to be easily cleaned.



Display and Control

Confidence and Control At Your Fingertips

Terumo System 1's Central Control Monitor allows complete control and monitoring of the entire perfusion case from a highly visible, central touch-screen display. Beneath the attractive, easy-to-view graphics lies an intelligently and intuitively designed software program for fast, simple configuration and running of a case.

- Device identification is fast and easy with the simple perfusion circuit graphic.
- Important patient parameters are always visible, even as menu screens are changed and viewed.
- Critical information and operating controls are accessible with one touch.
- All menu screens can be accessed with no more than two touches.
- Alarm, alert, status and error messages are easily viewable.
- Pulse, servo regulation to pressure or flow, master/follower operation for multi-ratio cardioplegia, cardioplegia volume tracking and cardioplegia dose delivery are easy to activate and control from the Central Control Monitor.

Creating your own perfusion display screen is quick and simple:

- Select a graphical orientation.
- Select and place the pumps on the graphical interface.
- Select and place the safety and monitoring devices on the graphical interface.
- Define responses to alert and alarm conditions.
- Assign displayed devices to physical components.
- Save the configuration.

Custom display screens are transferable between systems using a PC card. Access to perfusion display screens can be controlled by using the optional password security feature.



Intuitive Design

The display design uses a simple perfusion circuit graphic with icons for easy identification of system devices and parameters. Quickly access information or activate controls by touching the corresponding icon.




Device Controls and Display

Adjust pump speeds, the occluder position, FIO₂ setting, and gas flow rate using the control bars located on either side of the touch screen. View important patient parameters directly on the device icons.



Additional Functionality

Tabs along the bottom of the touch screen provide access to additional information and functionality without covering critical controls. Menu screens displayed at the top of the screen allow you to efficiently enter patient information, modify system settings, activate safety systems, monitor cardioplegia delivery, and view additional system status information.




Alert and Alarm Handling

Color-coded alert and alarm messages are displayed indicating the condition and status of the system. The device icon for the alert or alarm source flashes and the corresponding safety screen is automatically displayed so you can respond quickly.

Advanced Functionality with Reliability You Trust

Terumo System 1's pumps are smaller and offer more advanced features than the Sarns brand models that preceded them and continue to provide the same precision engineering and solid design that have made Sarns brand pumps legendary for reliability and dependability.




The roller pump raceway rotates in 15° increments, allowing you to minimize tubing lengths in your circuit.

Centrifugal Pumps

- Pumps can be controlled from the front panel knob or the Central Control Monitor.
- Remote-mounted pump motor helps optimize tube lengths.
- Non-invasive flow sensor eliminates need for a disposable flow probe.
- Integrates with system alerts and alarms.
- Advanced functions include:
 - Pulsatile flow operation (arterial pump only).
 - Servo regulation to maintain a constant flow or a constant positive or negative pressure.
 - User selectable "coast" response reduces the pump speed to 1500 RPM in response to an alert or alarm.



The centrifugal pump's compact size enables easy pole mounting of the entire system.



Self-adjusting tube clamps eliminate the need to replace tube inserts.

Roller Pumps

- Occlusion mechanism with audible feedback can be adjusted while the pump is running.
- Easy-to-use mounting hardware allows choice of base or pole mounting.
- Pumps can be controlled from the front panel knob or the Central Control Monitor.
- Front panel display shows pump operating status and any alerts or alarm messages affecting its operation.
- Preferred system direction for forward flow is user selectable.
- Advanced functions can be performed using any roller pump:

- Pulsatile flow operation (arterial pump only).
- Servo regulation to maintain a constant flow or a constant positive or negative pressure.
- Master/follower operation of two pumps for delivering multi-ratio cardioplegia.
- Cardioplegia volume or time doses can be set on the Central Control Monitor to stop the cardioplegia pump when delivery is complete.



The System 1 roller pumps are available in two sizes: 6" and 4" diameter.

Specifications

System 1 Base

	<u>801763</u>	<u>801764</u>
Voltage	100/115V 50/60 Hz	220/240V 50/60 Hz
Height	22.6 in (57.4 cm)	22.6 in (57.4 cm)
Width	35.2 in (89.4 cm)	35.2 in (89.4 cm)
Depth	26.5 in (67.3 cm)	26.5 in (67.3 cm)
Weight	262 lbs (118.8 kg)	262 lbs (118.8 kg)

Electronic O₂ Blender/Analyzer

Operating range:	<u>801188</u>
Flow	0 – 10 L/min
FiO ₂	0.21 – 1.00
Measured O ₂	21% – 100%

Central Control Monitor

	<u>802100</u>
Height	13 in (33 cm)
Width	15.7 in (39.9 cm)
Depth	3.0 in (7.6 cm)
Weight	15 lbs (6.8 kg)

Roller Pumps

	<u>801040 (small)</u>	<u>801041 (large)</u>
Pumphead diameter	4 in (10.2 cm)	6 in (15.2 cm)
Voltage	24VDC	24VDC
Height	12.5 in (31.8 cm)	12.5 in (31.8 cm)
Width	7.1 in (18.0 cm)	8.5 in (21.6 cm)
Depth	11.8 in (30.0 cm)	13.1 in (33.3 cm)
Weight	21 lbs (9.5 kg)	24 lbs (10.9 kg)
Operating range	0 – 4 L/min	0 – 10 L/min

Centrifugal Control Unit

	<u>801046</u>
Voltage	24VDC
Height	3.1 in (8.0 cm)
Width	7.3 in (18.4 cm)
Depth	8.5 in (21.6 cm)
Weight	2.4 lbs (1.1 kg)
Operating range	0 – 7 L/min

Flexible Halogen Lamps

	<u>801238</u>	<u>801558</u>
Size	33 in (83.8 cm)	15 in (38.1 cm)
Voltage	24VDC	24VDC

Functional Modules

	<u>Flow module 802018</u>	<u>Other modules</u>
Height	6.30 in (160.0 mm)	3.54 in (90.0 mm)
Width	3.63 in (92.2 mm)	1.06 in (27.0 mm)
Depth	1.33 in (33.8 mm)	3.03 in (77.0 mm)
Weight	0.95 lbs (0.43 kg)	0.27 lbs (0.12 kg)

Pressure 802112

Two pressure transducers per module
Operating range; (-250) mmHg – 900 mmHg
Maximum of 8 transducers

Temperature 802114

YSI 400, two temperature sensors per module
Operating range; 0 – 50° C
Maximum of 8 sensors

Flow 802018

Non-invasive flow measurement, one per module
Operating range; (-9.9) L/min to 9.9 L/min
Maximum of 4 modules

Ultrasonic 802111

Level Detection One alarm sensor, one alert sensor per module
Reservoir: Functions with hardshell reservoirs that have a wall thickness of 0.07 in – 0.15 in (1.8 mm – 3.8 mm)
Maximum of 1 module

Ultrasonic Air 802110

Bubble Detector One per module
Operating range;
3/8 in – 0.5 cc or larger up to 6 L/min
1/4 in – 0.3 cc or larger up to 3 L/min
Maximum of 4 modules

Electronic Venous 803480

Occluder One per module
Operating range; 0 – 100% flow on
1/4 in to 1/2 in tubing
Maximum of 1 module

Interface Module 802558

for CDI System 100 One CDI System 100 monitor per module
Maximum of 1 module (either CDI System 100 or CDI System 500)

Interface Module 803479

for CDI System 500 One CDI System 500 monitor per module
Maximum of 1 module (either CDI System 100 or CDI System 500)

Interface Module 802113

for RS-232 One CCPro Data system per module
Maximum of 1 module (either RS-232 or RS-485)

Interface Module 803518

for RS-485 One data management system per module
Maximum of 1 module (either RS-232 or RS-485)

Ordering Information

System 1 Base

1* 100/120V System 1	801763
or 220/240V System 1	801764
<i>* Shaded items are included with System 1</i>	
Programmed PC card system configuration	803739
2 Hand crank (includes 2 hand cranks)	801016
3 Hand crank bracket	802089
4 Central Control Monitor	802100

Base Options

5 Electronic O ₂ blender/analyzer*	801188
6 Pole-mounted blender*	164235

** Each blender requires a hose kit and adaptor set*

Hose kits

U.S. hose kit (3 hoses: green, yellow, black)	144194
Non-U.S. hose kit (3 hoses: blue, yellow, black)	144186

Hose adaptor sets

7 NCG hose adaptor set	144207
8 D.I.S.S. hose adaptor set	144215
9 Ohio Diamond hose adaptor set	144223

Center Poles

10 Crossbar fitting (required for each additional pole)	145980
11 2 ft (0.6 m) pole	16553301
3 ft (0.9 m) pole	131115
4 ft (1.2 m) pole	16553401

Flexible Halogen Lamp

12 33 in (83.8 cm) flexible halogen lamp	801238
15 in (38.1 cm) flexible halogen lamp	801558

Roller Pumps

13 Roller pump 6 in (15.2 cm) diameter	801041
14 Roller pump 4 in (10.2 cm) diameter	801040
15 Pole mount pump rest with bracket	801093

Integrated Centrifugal System

16 Drive motor	164267
17 Control unit	801046
Manual drive	164268
Pole mount centrifugal display bracket	804372
Centrifugal motor holder	7088

Sarns CCPro Data Perfusion Software

Data management software	194765
--------------------------------	--------

Flow Sensing

18 Flow module	802018
19 Non-invasive flow sensor 3/8 in (9.5 mm) ID x 3/32 in (2.4 mm) wall, reusable	6382
20 Mounting bracket (holds 2 modules)	801550

Level Detection

(One each included with System 1)

21 Level detect module	802111
22 Yellow transducer	195215
22 Red transducer	195274
Level sensor pads, gel included (60 per box)	195240

Air Bubble Detection

(One each included with System 1)

23 Air bubble detect module	802110
24 Cable assembly	149892
Pole clip sensor holder	149876

One sensor required per air bubble detection system

Ultrasonic air sensor (choose 1 of 3)

25 3/8 in x 3/32 in (9.5 mm x 2.4 mm)	5773
25 1/4 in x 3/32 in (6.4 mm x 2.4 mm)	5791
25 1/4 in x 1/16 in (6.4 mm x 1.6 mm)	5785
26 Air sensor bracket (optional)	5793

Venous Occluder

27 Occluder module	803480
28 Occluder head	806455

Pressure Monitoring

29 Pressure module	802112
30 Reusable pressure transducer	16433301
Pressure monitoring kit (10 per case)	16066100
Pressure transducer holder	22300030

Temperature Monitoring



31 Temperature module, YSI® series 400 compatible (2 temperatures per module)	802114
---	--------

YSI is a registered trademark of Yellow Springs Instruments

Data Interface Modules

32 Interface module for CDI System 100	802558
33 Interface module for CDI System 500	803479
34 Serial interface module RS-232	802113
35 Serial interface module RS-485	803518

Shaded items included with System 1




Terumo® Advanced Perfusion System 1

Sarns™
Module 5000
Console

Sarns™
7000 MDX
System

Sarns™
Modular
Perfusion System
8000

Sarns™
Perfusion System
9000



System 1

Terumo® Advanced
Perfusion System 1

Heir to the Sarns Legacy

Terumo System 1 may employ 21st century technologies, but it has inherited the best features from its Sarns ancestors, the longest commercially available line of heart-lung machines in the world. Since 1961, Sarns heart-lung machines have been legendary for their high-quality engineering, longevity, and customer satisfaction. Like those ancestors, Terumo System 1 is “built for life.”

Evolved Serviceability

System 1 is designed with a minimum of moving components, simplifying maintenance and repair. It has enhanced self-diagnostics. And its system modules can be replaced during use to minimize downtime if a problem should occur. A special software-based service tool was developed to speed diagnosis and reporting. All these features, supported by experienced technical and field service teams, help make System 1 the most confident choice in perfusion systems today.

Ask your Terumo Cardiovascular Systems representative for additional information about special service features and programs for System 1.



TERUMO CARDIOVASCULAR SYSTEMS CORPORATION
6200 Jackson Road, Ann Arbor, MI 48103-9300
USA
734 663 4145 phone
734 663 7981 fax
800 521 2818 toll free

TERUMO MEDICAL CORPORATION
Latin American Sales Division
6161 Blue Lagoon Drive, Suite 340, Miami, FL 33126
305 263 9835 phone
305 263 1219 fax

TERUMO EUROPE N.V.
Cardiovascular Division
Interleuvenlaan 40 B-3001, Leuven
Belgium
32 16 38 12 11 phone
32 16 40 02 49 fax

TERUMO CORPORATION
44-1, 2-chome, Hatagaya, Shibuya-ku, Tokyo 151-0072
Japan
81 3 3374 8111 phone
81 3 3374 8196 fax