Instructions for use - Mode d'emploi - Gebrauchsanleitung
Istruzioni per l'uso - Instrucciones de uso

INTEGRAL DRAINAGE SETS (IDS)

Ensembles de Drainage Intégré (IDS) - Integrale Drainage-Sets (IDS)
Set Integrali per Drenaggio (IDS) - Sets de Drenaje Integral (IDS)

MADE BY / Fabriqué par / Hergestellt von
Prodotto da / Manufacturado por
Integra NeuroSciences Implants S.A.
2905 Route des Dolines
06921 Sophia Antipolis cedex, **France**


Figure 1 - Figure 1 - Abb.1 - Figura 1 - Figura 1

Sterile - Sterilized with ethylene oxide gas. For single use only. Do not use open or damaged packages. Do not resterilize.

WARRANTY DISCLAIMER AND LIMITATION OF REMEDY

THERE IS NO EXPRESS OR IMPLIED WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR OTHER WARRANTY ON INTEGRA NEUROSCIENCES PRODUCTS DESCRIBED IN THESE INSTRUCTIONS FOR USE. UNDER NO CIRCUMSTANCES SHALL INTEGRA NEUROSCIENCES IMPLANTS BE LIABLE FOR MEDICAL EXPENSES OR ANY DIRECT, INDIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW. NO PERSON HAS ANY AUTHORITY TO BIND INTEGRA NEUROSCIENCES IMPLANTS TO ANY REPRESENTATION OR WARRANTY CONCERNING INTEGRA NEUROSCIENCES IMPLANTS PRODUCTS EXCEPT AS SPECIFICALLY SET FORTH HEREIN.

Descriptions and specifications appearing in Integra NeuroSciences printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Rx only.

Product returns to Integra NeuroSciences: Contact your local Representative for instructions on return of products.

DESCRIPTION

The Integral Drainage Set (IDS) consists of the following components (Figure 1- page 2):

1. A one-piece assembly, which includes:
 - Male Luer-lock connector (14) with self-sealing sampling port (13)
 - Small diameter drainage inlet tube (1) with attached clamp (12), check valve (2) and silicone band (17)
 - 700 ml graduated bag (8) [with permanent markings (7) and a drainage pressure reference arrow], 0.22µm hydrophobic vent filter (4), drip chamber (3), suspension string (5), with adjustable clip (6)
 - Large diameter drainage collector (9)
 - Outlet tube with clamp (10), watertight cap (11) and silicone band (17)
2. A 360° three-way stopcock (15) with two female Luer-lock connectors with cap (18) and 10 cm piece of tubing terminating in a male Luer-lock connector (16) with cap (19).

Note: Two three-way Stopcock Protection Boxes are supplied with some configurations.

Materials intended to contact body tissue and/or body fluids are silicone elastomer, silicone elastomer impregnated with barium sulfate, polyisoprene, polycarbonate, polyethylene, ethyl vinyl acetate (EVA), polyurethane and ruby.

INDICATIONS

The drainage sets are intended for cerebrospinal fluid (CSF) drainage, sampling and collection.

CONTRAINDICATIONS

The use of integral ventricular and lumbar drainage sets are contraindicated in patients with a history of blood-clotting abnormalities, or receiving anticoagulants, or with scalp or skin infections.

COMPLICATIONS

Possible complications include, but are not limited to: infection, obstruction, overdrainage, collapse of the ventricles, intracerebral hemorrhage, air leakage into the ventricle or subarachnoid space, and fluid leakage.

PRECAUTIONS

This product is recommended for single use only. *This product is sterilized with ethylene oxide. Do not use if the package is open or damaged. Use the device prior to the "Use Before" date on the package label.*

Caution - Do not resterilize. Integra NeuroSciences will not be liable for any direct, indirect, incidental or consequential damages resulting from or related to resterilization.

- Handle tubing with care. Before using forceps, cover the jaws with silicone elastomer tubing.
- The distal end of the lumbar catheter may be cut to desired length. Other modifications to the set are not recommended.

PREPARATION AND USE OF THE IDS [for Ventricular or Lumbar Drainage]:

Integra NeuroSciences recommends that the physician have duplicates of each sterilized item to be used during a procedure. Although these instructions are based on experience gathered to date, the physician may wish to vary procedure in accordance with clinical judgment.

Caution - Check all connections to ensure that fittings are tight and leak-free.

1. Using the suspension string and adjustable clip, hang the drainage bag to a support.
2. Position the reference arrow on the bag at the height appropriate for the desired drainage pressure (cm H₂O).
Note: To prevent backflow when the inlet tube clamp is open, the bag must be kept upright.
Note: To adjust drainage pressure, position the bag to the desired cm H₂O height. Use the scale on the bag to measure from the zero point (arrow on the bag) to the patient's eye level (ventricular drainage) or to the patient's spine (lumbar drainage).
Caution: To prevent unexpected bag movement, verify that the bag is attached securely once positioning is complete.
3. Occlude the drainage inlet tube with the clamp.
4. With the catheter in position, remove the cover from the male Luer-lock connector, leave the watertight cap in place. Attach the connector to the catheter adapter on the previously implanted catheter.
Caution: Overtightening may damage the connector.
Caution: Extra-care must be taken to isolate the assembly catheter/Luer-lock connector from any mechanical constraint.
Note: To disconnect the stopcock from the Y-connector, rotate the collar of the Y-connector counter-clockwise until the stopcock disengages.
5. To allow flow of CSF, release the clamps on the inlet tube and catheter.
Note: Regularly check system patency for any potential disconnection, obstruction ...
Note: Replace the IDS and individual components according to a hospital-approved schedule.

Protection from potential contamination:

Always proceed aseptically according to hospital procedures. To protect the three-way stopcock of the sampling port from potential contamination, use exclusively one of the following methods:

- Tight the cap on the Luer-lock connector. Observe aseptic precautions during sampling procedures.
- Use the Integra Stopcock Protection Boxes:
 - Soak the whole surface of the foams with either iodine polyvidone or chlorhexidine. Place the stopcock and close the box.
 - Add antiseptic every 8 hours or at each procedure.
 - Replace the box every 72 hours maximum.

Caution: Do not use other antiseptic solutions than polyvidone iodine or chlorhexidine, as they may damage the tubings.

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The wave logo and Integra NeuroSciences are trademarks of Integra LifeSciences Corporation.