Instructions for use  
**External Drainage Set (EDS)**

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**Sterile - Sterilized with ethylene oxide gas. For single use only. Do not use open or damaged packages. Do not resterilize.**

#### **WARRANTY DISCLAIMER AND LIMITATION OF REMEDY**

**INTEGRA NEUROSCIENCES IMPLANTS HEREBY DISCLAIMS ANY EXPRESS OR IMPLIED WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR OTHER WARRANTY ON INTEGRA NEUROSCIENCES IMPLANTS PRODUCTS DESCRIBED IN THESE INSTRUCTIONS FOR USE. UNDER NO CIRCUMSTANCES SHALL INTEGRA NEUROSCIENCES IMPLANTS BE LIABLE FOR MEDICAL EXPENSES OR ANY DIRECT, INDIRECT, INCIDENTAL, CONSEQUENTIAL, OR PUNITIVE DAMAGES. NO PERSON HAS ANY AUTHORITY TO MAKE ANY REPRESENTATION OR WARRANTY ON BEHALF OF INTEGRA NEUROSCIENCES IMPLANTS OTHER THAN THOSE SET FORTH HEREIN.**

**Descriptions and specifications appearing in Integra NeuroSciences Implants printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.**

**Caution:** Do not tip drainage cylinder while inside of cylinder is wet. Do not allow cylinder filter to become wet.  
**Rx only.**

Product returns to Integra NeuroSciences Implants: Contact your local Representative for instructions on return of Integra NeuroSciences Implants products.

## **DESCRIPTION**

The External Drainage Set (EDS) consists of the following components (Figure 1):

- An extension tubing line with, at one end a T fitting (male luer lock connector and latex-free access port), and at the other end a four-way stopcock (female luer lock connector).
- A one-piece drainage tube assembly, including a patient line with a T fitting (male luer lock connector and latex-free access port), a tubing (with clamps and four-way stopcock), an antireflux valve, a 75 ml cylinder (with hydrophobic filter, permanent marking for volume measurement and velcro straps for attachment) and a drainage tubing (with four-way stopcock latex-free access port and a male luer lock connector).
- A drainage bag with female luer lock connector, a suspension hook and a hydrophobic filter.
- A 350 mm H<sub>2</sub>O drainage pressure scale.

**Note** - Also available are Needle-Free EDS designed to be used with needleless syringes.

**Materials:** The product may contain one or more of the followings: acrylonitrile butadiene styrene (ABS), polyisoprene, silicone elastomer, polycarbonate, polyethylene, polyvinylchloride (PVC), and synthetic ruby.

## **INDICATIONS**

The drainage sets are intended for cerebrospinal fluid (CSF) drainage, sampling and collection.

## **CONTRAINDICATIONS**

External ventricular drainage is contraindicated for patients with a history of blood-clotting abnormalities, receiving anticoagulants or with scalp or skin infections.

## PRECAUTIONS

Do not use if the "Use by Date" on the package label has expired.  
Do not use open or damaged packages.  
Do not autoclave.  
Do not expose to temperatures exceeding 60°C (140°F ).  
Handle tubing with care. Before using forceps, cover the jaws with silicone elastomer tubing.

Modifications are not recommended.

## WARNINGS

These devices are intended for one use only. All components are extremely difficult to clean after being exposed to biological materials and adverse patient reactions may result from reuse. **Integra NeuroSciences Implants will not be liable for any direct, indirect, incidental consequential, or punitive damages resulting from reuse of the product.**

## COMPLICATIONS

Possible complications include but are not limited to infection, obstruction, overdrainage, collapse of the ventricles, intracerebral hemorrhage, air leakage into the ventricle or subarachnoid space, and fluid leakage.

## STERILIZATION

These devices are sterilized with ethylene oxide gas and are supplied in sealed packages designed to protect the product during normal shipment and handling.

**Caution - Do not resterilize. Integra NeuroSciences Implants will not be liable for any direct, indirect, incidental, consequential or punitive damages resulting from or related to resterilization.**

## REPLACEMENT CREDIT

To assure that the packaged product is received in the same condition in which it left the manufacturer, the container is sealed after final inspection. Replacement credit cannot be allowed if the seal is broken, or if the "Use by Date" has expired.

## PREPARATION AND USE OF THE EDS [for Ventricular or Lumbar Drainage]:

It is recommended that the physician have duplicates of each sterilized item to be used during a procedure. Thus, if aseptic technique is compromised, the procedure can continue.

Although these instructions are based on experience gathered to date, the physician may wish to vary procedure in accordance with clinical judgment.

**Note** - The "EDS Easy Stand" from Integra NeuroSciences Implants may be used to support the EDS and allows for its quick set up. Stopcock Protection Boxes are also available separately.

**See figures 2 and 3:**

1. Proceed with the installation of the EDS Easy Stand following its Instructions for Use, or attach the 350 mm drainage pressure scale, with self-adhesive backing, to a support (such as an IV pole). Place the zero point of the scale at the patient's eye level (ventricular drainage) or spine level (lumbar drainage).

**Caution** - If the level of the patient's head (ventricular drainage) or spine (lumbar drainage) is changed, readjust the scale to the new zero point.

2. Insert the male connector at the end of the drainage tube on the cylinder assembly into the female connector on the drainage bag.
3. Hang the drainage bag from a suitable support well below the level of the patient.
4. Position the cylinder against the support with the reference arrow at the desired drainage pressure height. Pull the cylinder velcro straps tightly around the support and firmly press the strap loop against the cylinder. Firmly pull on the cylinder (up and down) to verify the integrity of its positioning.

When using the EDS Easy Stand, snap the EDS cylinder into the Easy Stand cylinder holder. Adjust the cylinder to the desired drainage pressure height and secure it in position, using the thumb screw. Align the reference arrow on the cylinder with the desired drainage pressure. Secure the cylinder to the cylinder holder with one velcro strap.

**Caution** - Verify the safe positioning of the cylinder.

**Caution** - To prevent backflow, the cylinder and the check valve must be kept upright, and the tubing must be secured in the velcro straps.

**Note** - To change the drainage pressure, reposition the cylinder to the new desired mm H<sub>2</sub>O height. Use the scale to measure from the zero point up to the arrow on the cylinder.

5. Occlude the patient line with the clamp.
6. Attach the Luer lock T-connector on the patient line extension tubing to the catheter adapter on the previously implanted catheter.

**Note** - To disconnect the catheter adapter from the T-connector, rotate the collar of the T-connector counterclockwise until the catheter adapter disengages.

**Caution** - Before draining, check all connections to ensure that fittings are tight and leak-free. Do not overstress connections. All connections on the one-piece drain tube assembly are permanently glued and cannot be disconnected without breakage.

7. Rotate the four-way stopcocks on the extension tubing and patient line to their drainage positions. To allow flow of CSF, release the clamp on the inlet tube. Verify adequate CSF flow.

8. Select the desired four-way stopcock position under the cylinder for (see Figure 4):
  - collection of CSF in the graduated cylinder
  - collection of CSF in the graduated cylinder and sampling of CSF through the lower sampling port
  - drainage of CSF into the drainage bag
  - drainage of CSF into the drainage bag and sampling of CSF through the lower sampling port

**Caution:**

- Simultaneous drainage and pressure monitoring (stopcock in "all opened" position) is not recommended for precise pressure measurements.
- Do not overfill the cylinder. Keep the cylinder and the check valve upright at all times.
- To empty the CSF collection in the drainage bag, first occlude the patient line with the clamp then drain the cylinder. Release the clamp once drainage to the bag is finished.
- Patient line should be occluded with the clamp during any change of patient's position. Reset the zero level prior to unclamp the line.

**Notes** - Regularly check system patency for any potential disconnection, obstruction, or other problem. Replace the EDS and individual components according to a hospital-approved schedule.

**Use of the Needle-Free EDS sampling port:**

The port is a two-way valve, completely closed system, activated by the Luer taper of a syringe.

1. Remove and throw away the cap. Be careful not to touch the opening of the valve.
2. Screw the needleless syringe into the valve. The valve is now open.
3. Proceed with sampling.
4. Remove syringe. The valve closes when the syringe is withdrawn. **Place a new sterile cap on the valve after each procedure.**

Always proceed aseptically according to hospital procedures.

**Protection from potential contamination:**

Always proceed aseptically according to hospital procedures.

To protect the four-way stopcock of the access port from potential contamination, use exclusively one of the following methods:

- Tighten the cap on the Luer-lock connector. Observe aseptic precautions during sampling procedures.
- Use the Integra Stopcock Protection Boxes:
  1. Soak the whole surface of the foams with antiseptics (iodine polyvidone, chlorhexidine, etc.). Place the stopcock and close the box.
  2. Add antiseptic every 8 hours or at each procedure.
  3. Replace the box every 72 hours maximum.

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The wave logo and Integra NeuroSciences are trademarks of Integra LifeSciences Corporation.