

Expanding the Treatment Options for Atrial Fibrillation

Temperature sensor provides continuous monitoring of cooling effectiveness

EXTEND THE POTENTIAL OF IRRIGATED

ABLATION WITH ROBOTIC CATHETER CONTROL

- Integrated irrigated ablation catheter compatible with the Sensei® X Robotic Catheter System
- Leverages the navigation capability of the Artisan Extend™ Control Catheter in a smaller, more flexible profile

THE BENEFITS YOU EXPECT FROM IRRIGATED CATHETER TECHNOLOGY

- Designed to reduce char and to deliver larger, deeper lesions than non-irrigated ablation catheters
- Developed to increase the effectiveness of complex cardiac arrhythmia treatment, including atrial fibrillation
- Compatible with the Stockert™ 70 or Stockert™ EP Shuttle RF Generators and the CoolFlow™ 70 Pump

Flexible low-profile sheath and ablation catheter with six degrees of articulation facilitate access to difficult-to-reach anatomy



LYNX CATHETER LESION DATA

Contact Force 20g*

Power settings (W)	15	20	25
Average power output (W)	13	18	23
Maximum temperature (°C)	38	39	39
Lesions formed		£ .	243

Contact Force 30g*

Power settings (W)	15	20	25
Average power output (W)	13	18	22
Maximum temperature (°C)	39	40	41
Lesions formed		•	

^{*} Lesion examples in vitro with chicken breast meat in a wet lab setting. The exact lesion results cannot be translated directly to the in-vivo experience. These lesions were formed using the power control mode for 60 seconds each.

SMALL. FLEXIBLE. IRRIGATED

The Lynx™ Robotic Ablation Catheter integrates a robotically steerable irrigated ablation catheter into a robotically steerable sheath to provide a high level of control, flexibility, and stability during lesion formation.



Ablation catheter articulation 0–270° in any direction with a variable bend radius

TIP ELECTRODE	NUMBER OF ELECTRODES		ELECTRODE DIAMETER	TOTAL EFFECTIVE LENGTH	TOTAL ARTICULATING LENGTH
4 mm	2	2 mm	7F	82 cm	10 cm



www.hansenmedical.com

Corporate Headquarters

800 East Middlefield Road Mountain View, CA 94043 Telephone: 650 404 5800 Toll Free: 888 404 5801 Fax: 650 404 5901

Hansen Medical EMEA

Boston House Business Centre 69-75 Boston Manor Road, Brentford Middlesex, TW8 9JJ, United Kingdom Telephone: +44 208 847 8608 Fax: +44 208 847 8669

This device is not for sale in the United States.

The Lynx™ Robotic Ablation Catheter requires CE mark approval and is not commercially available in the EU.