# FemoSeal®

Vascular Closure System

ENGLISH  
SVENSKA  
DEUTSCH  
FRANÇAIS  
DANSK  
SUOMI  
NEDERLANDS  
NORSK  
ITALIANO  
ESPAÑOL  
PORTUGUESE  
MAGYAR  
ČEŠTINA  
POLSKI  
ΕΛΛΗΝΙΚΑ  
EESTI  
TÜRKIYE'YE İLİŞKİN

REF 11202

**RADI**  
Improving Patient Care


# CONTENTS

ENGLISH	3
SVENSKA	9
DEUTSCH	15
FRANÇAIS	22
DANSK	29
SUOMI	35
NEDERLANDS	41
NORSK	48
ITALIANO	54
ESPAÑOL	61
PORTUGUESE	68
MAGYAR	75
ČEŠTINA	82
POLSKI	88
ΕΛΛΗΝΙΚΑ	95
EESTI	102
TÜRKIYE'YE İLIŞKİN	108

# ENGLISH

## Device description

FemoSeal® Vascular Closure System is a medical device indicated for closure of an arterial puncture after percutaneous catheterization through the femoral artery. FemoSeal® Closure Elements consist of two resorbable polymer discs, the Inner Seal and Outer Locking Disc, which are held together by a resorbable multifilament, thereby mechanically sealing the arteriotomy. Hemostasis is achieved by mechanical means.

## Indication for use

FemoSeal® Vascular Closure System is indicated for use in closing the femoral arterial puncture (arteriotomy) in patients who have undergone percutaneous catheterization using a 7F (2.33 mm) or smaller procedural sheath.

## Contraindications

FemoSeal® Vascular Closure System is contraindicated in patients with arteriotomies in which sheaths or devices larger than 7F (2.33 mm) have been used.

## Warnings

- Do not use if the package has been damaged or any sterile barrier is not intact.
- Do not use after expiry date - the biodegradable components may not perform adequately.
- Do not use if any items appear damaged or defective in any way.
- Use of FemoSeal® Vascular Closure System where bacterial contamination of the procedural sheath or the surrounding tissue may have occurred, may cause infection.
- If it is suspected that the posterior arterial wall has been punctured or more than one arterial puncture has been made, do not rely solely on FemoSeal® Vascular Closure System to achieve arterial hemostasis. Use additional manual or mechanical compression.
- If the puncture site is at, or distal to the bifurcation of the femoral artery,

FemoSeal® Vascular Closure System should not be used due to the risk of the Inner Seal being positioned incorrectly. This event may result in bleeding complications and/or disruption to normal blood flow.

- If there is persistent arterial bleeding from the incision site, do not cut the multifilament until hemostasis is achieved. If hemostasis is not achieved, fasten the multifilament with a sterile wound dressing and apply supplementary compression until hemostasis is achieved. In the case of persistent arterial bleeding, significant bleeding complications may occur which could result in patient injury or death.

### **Special patient populations**

The safety and effectiveness of FemoSeal® Vascular Closure System has not been established in the following patient populations:

- Patients with pre-existing autoimmune disease.
- Patients undergoing therapeutic thrombolysis.
- Patients with clinically significant peripheral vascular disease.
- Patients with uncontrolled hypertension (> 220 mmHg systolic or > 110 mmHg diastolic).
- Patients with a bleeding disorder, including thrombocytopenia (< 100,000 platelet count), or anemia (Hgb < 10 mg/dl).
- Patients having an inner lumen of the common femoral artery smaller than 5 mm.
- Patients with a vascular graft or stent at the puncture site.
- Patients who are pregnant or lactating.
- Patients below the age of 18 years.

### **Precautions**

- FemoSeal® Vascular Closure System deployment procedure should be performed by physicians with adequate training in the use of the device.
- Perform a limited femoral angiogram prior to FemoSeal® Vascular Closure System deployment.
- Discontinue procedure if:
  - Lumen diameter of common femoral artery < 5 mm.
  - Stenosis and/or significant plaque present in the vicinity of the femoral arterial puncture site.
  - Arterial puncture is at, or distal to, the common femoral artery bifurcation.
  - Anomalous branches or vessel abnormalities present in the vicinity of the femoral arterial puncture site.
- If any hematoma is present, extra care must be taken for correct insertion of the sheath into the artery.
- Use a single wall puncture technique. Do not puncture the posterior wall of the artery.
- Observe sterile technique at all times when using FemoSeal® Vascular Closure System.
- FemoSeal® Vascular Closure System is for single use only and should not be

- resterilized or reused in any manner.
- If deployment of the Inner Seal meets unexpected resistance, discontinue the procedure.
- For correct deployment of Inner Seal, reposition your thumb so that the button can spring back freely.
- Ensure the tip of the sheath of FemoSeal® Unit is under the skin surface before deployment of the Outer Locking Disc. The Outer Locking Disc may inadvertently be deployed above skin level in patients with a short distance between the femoral artery and the skin level.
- For correct deployment, a skin incision may be necessary before deployment.
- If the Inner Seal is inside the artery, but it is not possible to fully depress the button and deploy the Outer Locking Disc, surgery may be required to remove FemoSeal® Closure Elements.
- If repuncture of the same femoral artery becomes necessary within 90 days, repuncture should be made at least one centimeter proximal to the previous FemoSeal® Vascular Closure System access site.
- Instruct the patient to follow physician's orders regarding closure site inspection.
- Instruct the patient to carry the Patient Information Card for the next 90 days.

### **Precautions at time of discharge**

Before considering discharge, assess the patient for the following clinical conditions:

- Bleeding and/or hematoma at the closure site.
- Pain while walking.
- Signs of infection at the closure site.

### **Adverse events**

Adverse events that may occur and suggested treatment:

- *Bleeding and hematoma* – Pressure may be applied to the puncture site using digital or manual pressure, or using a compression device such as FemoStop® Femoral Compression System (Radi Medical Systems) to provide supplementary compression.
- *AV fistula or pseudoaneurysm* – If suspected, the condition may be evaluated with duplex ultrasound. When indicated, ultrasound guided compression of a pseudoaneurysm, e.g. FemoStop® Femoral Compression System (Radi Medical Systems), may be used. For severe bleeding complications, other therapeutic treatment may be appropriate.
- *FemoSeal® Closure Element non-deployment* – If FemoSeal® Closure Elements pull out with FemoSeal® Unit upon withdrawal, apply manual or mechanical pressure per standard procedure. Examine the device to ensure all resorbable components have been withdrawn.
- *Ischemic symptoms, arterial obstruction* – Should ischemic symptoms appear, consider ultrasound examination to locate position of Inner Seal since

vascular surgery may be indicated. Consider thrombolytic treatment. Monitor the patient for at least 24 hours.

- *Infection* – Any sign of infection at the puncture site should be taken seriously and the patient monitored carefully. Consider antibiotic treatment.
- *Edema* - Resorption of the implanted components is associated with a slight foreign body reaction. On rare occasions this can lead to a local edema.

## **FemoSeal® Vascular Closure System insertion procedure**

The procedure is described in the following steps:

### **1. Assess puncture site location by fluoroscopy.**

Perform a limited femoral angiogram prior to FemoSeal® Vascular Closure System deployment.

**CAUTION:** Discontinue procedure if:

- Lumen diameter of common femoral artery < 5 mm.
- Stenosis and/or significant plaque present in the vicinity of the femoral arterial puncture site.
- Arterial puncture is at, or distal to, the common femoral artery bifurcation.
- Anomalous branches or vessel abnormalities present in the vicinity of the femoral arterial puncture site.

**CAUTION:** If any hematoma is present, extra care must be taken for correct insertion of the sheath into the artery.

### **2. Position the sheath of FemoSeal® Unit inside the artery.**

- a) Insert FemoSeal® Dilator into FemoSeal® Unit until they snap together and marked signs meet (see figures 1a, 1b).

**NOTE:** Grip FemoSeal® Dilator near its tip, advance carefully in small increments.

- b) Insert the guidewire through the procedural sheath.

- c) Remove the procedural sheath and hold the guidewire in position (see figure 2).

- d) Thread FemoSeal® Unit over the guidewire and insert the sheath into the artery (see figures 3a, 3b).

**NOTE:** Leave approximately 25 cm of the guidewire visible.

### **3. Grip and squeeze the wings of FemoSeal® Safety Catch, lift and retract FemoSeal® Safety Catch slowly along with the guidewire until completely removed (see figures 4a, 4b, 4c, and 4d).**

**NOTE:** Hold the device in position with the opposite hand. When FemoSeal®

Safety Catch with dilator has been retracted a few centimeters (2 - 4 cm), blood should appear in the proximal part of the sheath, verifying that the sheath is in the artery.

**NOTE:** If no blood appears, retraction may be disrupted. Reintroduce the dilator and then the sheath into the artery and repeat the steps from 3.

- 4. Deploy the Inner Seal by pressing down the button completely in one continuous movement (see figures 5a, 5b).**

**CAUTION:** If deployment of the Inner Seal meets unexpected resistance, discontinue the procedure.

- 5. Pull back FemoSeal® Unit until the depressed button springs back (see figure 6).**

**CAUTION:** For correct deployment of Inner Seal, reposition your thumb so that the button can spring back freely.

**NOTE:** The button springs back once the Inner Seal makes contact with the inside of the arterial wall and a resistance is generated.

**NOTE:** The Inner Seal brings about the hemostasis.

- 6. Deploy the Outer Locking Disc by completely pressing down the button once again (see figure 7).**

**CAUTION:** Ensure that the tip of the sheath of FemoSeal® Unit is under the skin surface before deployment of the Outer Locking Disc. The Outer Locking Disc may inadvertently be deployed above skin level in patients with a short distance between the femoral artery and the skin level. For correct deployment, a skin incision may be necessary before deployment.

**NOTE:** A steep puncture angle shortens the distance between the femoral artery and the skin level.

**NOTE:** Keep a slight tension in the multifilament during deployment of the Outer Locking Disc until the button is completely pressed down. This step locks the Inner Seal against the arterial wall.

**CAUTION:** If the Inner Seal is inside the artery, but it is not possible to fully depress the button and deploy the Outer Locking Disc, surgery may be required to remove FemoSeal® Closure Elements.

- 7. Remove FemoSeal® Unit (see figure 8).**

## 8. Verify hemostasis.

**WARNING:** If there is persistent arterial bleeding from the incision site, do not cut the multifilament until hemostasis is achieved. If hemostasis is not achieved, fasten the multifilament with a sterile wound dressing and apply supplementary compression until hemostasis is achieved. In the case of persistent arterial bleeding, significant bleeding complications may occur which could result in patient injury or death.

## 9. Cut the multifilament below skin level using a sterile instrument (see figure 9).

### Storage, packaging and disposal

- Sterile and non-pyrogenic in unopened and undamaged package.
- FemoSeal® Vascular Closure System should be stored in a cold location (room temperature not exceeding 25 °C).
- FemoSeal® Vascular Closure System contains resorbable materials that degrade by exposure to heat and moisture; therefore, the device may not be resterilized.
- Dispose of contaminated units, components, and packaging materials utilizing standard hospital procedures and universal precautions for biohazardous waste.

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Radi Medical Systems AB