


St. Jude Medical Infinity™ Directional Lead

Models: 6170 | 6171 | 6172 | 6173 | 6178 | 6179 | 6180 | 6181



Compatibility

St. Jude Medical Infinity™ IPG (Model 6660)

St. Jude Medical Infinity™ IPG (Model 6662)

Lead Specifications

Model #	Lead			Electrode		
	Length	Diameter	Marker Band	Length	Longitudinal Spacing	Configuration
6170	30 cm	1.27 mm	Black	1.5 mm	0.5 mm	1-3-3-1
6172	40 cm	1.27 mm	Black	1.5 mm	0.5 mm	1-3-3-1
6178	30 cm	1.27 mm	White	1.5 mm	0.5 mm	1-3-3-1
6180	40 cm	1.27 mm	White	1.5 mm	0.5 mm	1-3-3-1
6171	30 cm	1.27 mm	Black	1.5 mm	1.5 mm	1-3-3-1
6173	40 cm	1.27 mm	Black	1.5 mm	1.5 mm	1-3-3-1
6179	30 cm	1.27 mm	White	1.5 mm	1.5 mm	1-3-3-1
6181	40 cm	1.27 mm	White	1.5 mm	1.5 mm	1-3-3-1

DBS Extensions with eXtend™ Lead Technology

Model #	Length	Marker Band
6371	50 cm	Black
6372	60 cm	Black
6373	90 cm	Black
6377	50 cm	White
6378	60 cm	White
6379	90 cm	White

Accessories

Model #	Description
6010	Guardian™ Burr Hole Cover
3014	Multilead Trial Cable



Rx Only

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

CE marking does not necessarily indicate regulatory approval status for all markets. Please refer to the instructions for use for a full listing of indications, contraindications, warnings and precautions.

Unless otherwise noted, TM indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL and the nine-squares symbol are trademarks and service marks of St. Jude Medical, Inc. and its related companies. © 2015 St. Jude Medical, Inc. All Rights Reserved.

St. Jude Medical Inc.

Global Headquarters
One St. Jude Medical Drive
St. Paul, MN 55117
USA
T +1 651 756 2000 | **F** +1 651 756 3301

St. Jude Medical S.C., Inc.

Americas Division
6300 Bee Cave Road
Bldg. Two, Suite 100
Austin, TX 78746
USA
T +1 512 286 4000 | **F** +1 512 732 2418

SJM Coordination Center BVBA

The Corporate Village
Da Vincilaan 11-Box F1
B-1935 Zaventem, Belgium
T +32 2 774 68 11 | **F** +32 2 772 83 84

St. Jude Medical Brasil Ltda.

Rua Itapeva, 538
5º ao 8º andares
01332-000 – São Paulo – SP
Brazil
T +55 11 5080 5400 | **F** +55 11 5080 5423

St. Jude Medical (Hong Kong) Ltd.

Suite 1608, 16/F Exchange Tower
33 Wang Chiu Road
Kowloon Bay, Kowloon
Hong Kong SAR
T +852 2996 7688 | **F** +852 2956 0622

St. Jude Medical Japan Co., Ltd.

Shiodome City Center 15F
1-5-2 Higashi Shinbashi, Minato-ku
Tokyo 105-7115
Japan
T +81 3 6255 6370 | **F** +81 3 6255 6371

St. Jude Medical Australia Pty, Ltd.

17 Orion Road
Lane Cove, NSW 2066
Australia
T +61 2 9936 1200 | **F** +61 2 9936 1222

SJMprofessional.com