

CONNECTION CABLES

MICRODEEP® DEPTH AND FORAMEN OVALE ELECTRODES



CONNECTION CABLES CBL**

- ◆ 5 to 18 contacts
- ◆ Length of cable available : 30 to 180 cm
- ◆ Only one electrode can be connected to one cable
- ◆ Touch-proof connectors (DIN 42802)

Microdeep® electrodes	Ref.	*Length of cable
D08-05AM	CBL-05/*	
D08-08AM	CBL-08/*	
D08-10AM	CBL-10/*	
D08-12AM	CBL-12/*	* 30 cm
D08-15AM	CBL-15/*	* 120 cm
D08-15BM	CBL-15/*	* 150 cm
D08-15CM	CBL-15/*	* 180 cm
D08-18AM	CBL-18/*	
D08-18CM	CBL-18/*	

Foramen ovale electrodes	Ref.	*Length of cable
F08-05AM-K #	CBL-05/L	* 30 cm * 120 cm
F08-10AM-K #	CBL-10/L	* 150 cm * 180 cm

Bold character references correspond to most current models available in stock.

EXTRA GREEN CONNECTION HOLDER SCC-F08S-10

- ◆ Spare part for the foramen electrode **F08****
- ◆ Single-use sterilised item, individually packed, delivered per lot of 10

