

GEISTER



The better way to operate™



New

CARDIOFRAME™ DR

Now with a lower profile and deeper reach than ever!

Universal sternum retractor system and operating platform for valve surgery and beating heart surgery.

CARDIOFRAME™ DR

This low profile operating platform keeps operating area free for surgeon's maneuvers and combines maximum comfort through a unique quick snap-in coupling for accessories with less traumatic tissue handling, stable fixation and excellent exposure.

New malleable hooks provide for a deeper reach as well as for versatile adjustment to individual requirements and patient anatomy. Carriers can be moved and fixed to any desired position on arms or spreading bar.

The wide range of GEISTER valve and beating heart accessories can also be used with this platform. The system also accommodates some of the most common disposable arms for beating heart surgery.

Quick change of blades

Order information:

- 29-0955 CardioFrame™ Retractor, Frame (29-0955.BD) with 2 inserts (29-0955.IN) and 2x Cooley blades 30mm deep (29-0956.30)
- 29-0940 Sliding Carrier, snap-in coupling (3 recommended)
- 29-0955.30 4x Morse blades 30mm deep
- 29-0955.40 4x Morse blades 40mm deep
- 29-0956.30 2x Cooley blades 30mm deep (included in 29-0955)
- 29-0956.40 2x Cooley blades 40mm deep

29-0940

- 29-1430.CF Cooley Atrial Hook small 20mm, malleable shaft
- 29-1431.CF Cooley Atrial Hook medium 35mm, malleable shaft
- 29-1432.CF Cooley Atrial Hook large 45mm, malleable shaft

29-0955.40

- 29-0941 Valve Hook small 18x52mm, short 14cm shaft
- 29-0942 Valve Hook medium 18x80mm, short 14cm shaft
- 29-0943 Valve Hook large 22x80mm, short 14cm shaft

29-0956.40

- 29-0941L Valve Hook small 18x52mm, long 19cm shaft
- 29-0942L Valve Hook medium 18x80mm, long 19cm shaft
- 29-0943L Valve Hook large 22x80mm, long 19cm shaft

- 29-0944L Valve Hook small 18x52mm, long 19cm curved shaft
- 29-0945L Valve Hook medium 18x80mm, long 19cm curved shaft
- 29-0946L Valve Hook large 22x80mm, long 19cm curved shaft

GEISTER®

The better way to operate™

Geister Medizintechnik GmbH
 Föhrenstraße 2
 D-78532 Tuttlingen/Germany
 Phone +49-(0)74 61-96 62 40
 Telefax +49-(0)74 61-96 62 422
 E-mail: info@geister.com
 Internet: www.geister.com

Distributor