

# ACross™ Transseptal Access System

## 8.5 F

### Product Highlights

- Directional fins orient each component's curve in the same plane
- Needle management system with interlocking components controls the needle deployment depth
- Braided sheath with ultra soft tip
- 8.5 F sheath ID accommodates a wide variety of catheters
- Two curve options available (LAMP™ 90, SL1™)

### Ordering Information


Contents: 8.5 F Swartz™ Braided Transseptal Guiding Introducer, Transseptal Dilator, Transseptal Needle, Guidewire and Stylet (1 unit per box)

ACross is packaged in a single tray that contains 4 separate components.  
(1. Guidewire, 2. Sheath, 3. Dilator, 4. Needle and Stylet)

Reorder Number	Sheath French Size	Sheath Usable Length (cm)	Sheath Curve Type	Needle Curve Type
406977	8.5	63	LAMP™ 90	BRK™
406978	8.5	63	SL1™	BRK™
406980	8.5	81	SL1™	BRK™
406981	8.5	63	LAMP™ 90	BRK-1™

#### Separate components:

Figure 1.



#### Assembly shown here:

Figure 2.

