

CPS Direct™ PL

Peelable Outer Guide Catheter

Product Highlights

- Unique SiteMark™ tungsten marker stripes provide superior fluoroscopic visibility to verify torque transfer
- Compatible with CPS Aim™ inner catheter and CPS Luminary™ bideflectable catheter with lumen to modify shape and extend reach if necessary
- EvenPeel™ stripes provide more smooth and reliable peeling for worry-free sheath removal



Ordering Information

Included: sheath with hemostasis valve attached, dilator and 2 valve bypass tools

Model Number	Curve Shape	Available Length (cm)	Overall Length (cm)	Inner Diameter (F/mm)	Outer Diameter (F/mm)
410210	Straight (OC-STR)	47	50,7	7/2,44	9/3,00
410211	Multipurpose (OC-MP)	47	50,7	7/2,44	9/3,00
410212	115° (OC-115)	47	50,7	7/2,44	9/3,00
410213	135° (OC-135)	47	50,7	7/2,44	9/3,00
410214	Wide (OC-W)	47	50,7	7/2,44	9/3,00
410215	Extra Wide (OC-XW)	47	50,7	7/2,44	9/3,00
410216	Right Sided (OC-R)	47	50,7	7/2,44	9/3,00
410224	145° (OC-145)	47	50,7	7/2,44	9/3,00
410217	Straight (OC-STR)	54	57,7	7/2,44	9/3,00
410218	Multipurpose (OC-MP)	54	57,7	7/2,44	9/3,00
410219	115° (OC-115)	54	57,7	7/2,44	9/3,00
410220	135° (OC-135)	54	57,7	7/2,44	9/3,00
410221	Wide (OC-W)	54	57,7	7/2,44	9/3,00
410222	Extra Wide (OC-XW)	54	57,7	7/2,44	9/3,00
410223	Right Sided (OC-R)	54	57,7	7/2,44	9/3,00
410225	145° (OC-145)	54	57,7	7/2,44	9/3,00



Separately Available Accessories

Model Number	Name	Type
410194	CPS Direct™ PL Valve Bypass Tool (Pack of 2)	Valve bypass tool
410195	CPS Direct™ PL Inner Catheter SafeSheath™ Sealing Adapter	Self-sealing valve
410190	CPS™ Implant Kit (Includes Needle, Syringe and 0,035" Guidewire)	Implant Kit

Customer Support: 46-8-474-4756

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use. Devices depicted may not be available in all countries. Check with your St. Jude Medical representative for product availability in your country. Unless otherwise noted, ™ indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2011 St. Jude Medical, Inc. All rights reserved.

CPS Direct™ SL II

Slittable Outer Guide Catheter

Product Highlights


- Integrated hub and hemostasis valve
- Increased curve retention and optimized catheter body structure for improved kink resistance
- Soft tip to lessen risk of traumatic insertion



Ordering Information

Included: dilator and 2 valve bypass tools

Model Number	Curve Shape	Available Length (cm)	Overall Length (cm)	Inner Diameter (F/mm)	Outer Diameter (F/mm)
DS2C001	Straight	47	50,7	7/2,44	9/3,00
DS2C002	115°	47	50,7	7/2,44	9/3,00
DS2C003	135°	47	50,7	7/2,44	9/3,00
DS2C004	Wide	47	50,7	7/2,44	9/3,00
DS2C005	X-Wide	47	50,7	7/2,44	9/3,00
DS2C006	Right Side	47	50,7	7/2,44	9/3,00
DS2C011	Straight	54	57,7	7/2,44	9/3,00
DS2C012	115°	54	57,7	7/2,44	9/3,00
DS2C013	135°	54	57,7	7/2,44	9/3,00
DS2C014	Wide	54	57,7	7/2,44	9/3,00
DS2C015	X-Wide	54	57,7	7/2,44	9/3,00



Separately Available Accessories

Model Number	Name	Type
DS2A003	CPS™ Universal Slitter	Slitter
DS2A004	CPS Direct™ SL Valve Bypass Tool	Valve bypass tool

Customer Support: 46-8-474-4756

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use. Devices depicted may not be available in all countries. Check with your St. Jude Medical representative for product availability in your country. Unless otherwise noted, ™ indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2011 St. Jude Medical, Inc. All rights reserved.

CPS Aim™ SL

Slittable Inner Catheter Cannulator with Integrated Valve


Product Highlights

- Integrated hemostasis valve in slittable catheter design
- Hydrophilic coating on outer catheter surface to enable smooth passage
- New catheter shaft/braid pattern for a kink-resistant and torqueable cannulator



Ordering Information

Model Number	Curve Shape	Available Length (cm)	Overall Length (cm)	Inner Diameter (F/mm)	Outer Diameter (F/mm)
DS2N024	CSL	65	68	5/1,83	7/2,29
DS2N025	AL2	65	68	5/1,83	7/2,29



Separately Available Accessories

Model Number	Name	Type
DS2A001	CPS Aim™ SL Inner Catheter SafeSheath™ Sealing Adapter	Self-sealing valve
DS2A002	CPS Aim™ SL Valve Bypass Tool	Valve bypass tool
DS2A003	CPS™ Universal Slitter	Slitter

Customer Support: 46-8-474-4756

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use. Devices depicted may not be available in all countries. Check with your St. Jude Medical representative for product availability in your country. Unless otherwise noted, ™ indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2011 St. Jude Medical, Inc. All rights reserved.

CPS Courier™

Guidewires

Product Highlights

- Helps physicians more easily subselect the target coronary branch vein and deliver the LV lead to its preferred destination



Ordering Information


Model Number	Distal Support	Length (cm)	Units per box	Diameter (in)
DS2G001	Soft	195	5	0,014
DS2G002	Medium	195	5	0,014
DS2G003	Firm	195	5	0,014
DS2G004	Extra Firm	195	5	0,014

CPS Duo™

Stylet Guidewire System

Product Highlights

- Enables optimal subselection of the branch vein and offers greater maneuverability and control of the LV lead



Ordering Information

Model Number	Type	Lengths (cm)	Diameter
DS2M001	CPS Duo™ Stylet	75; 86	OD: 0,014" LV lead lumen compatible ID: 0,012" compatible
DS2M006	CPS Duo™ Guidewire	195	0,012"

Customer Support: 46-8-474-4756

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use. Devices depicted may not be available in all countries. Check with your St. Jude Medical representative for product availability in your country. Unless otherwise noted, ™ indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2011 St. Jude Medical, Inc. All rights reserved.

CPS Luminary™

Bideflectable Catheter with Lumen



Product Highlights

- Two deflectable curves - a large deflection to facilitate coronary sinus cannulation and a small deflection for target vein subselection
- Soft, atraumatic tip features bipolar mapping electrodes to confirm CS entry
- CPS Direct™ SL outer guide catheters can be tracked over CPS Luminary™ catheter into the target vein

Ordering Information

Model Number	Description	Working Size (cm)	Overall Length (cm)	Size (F)	Inner Lumen Diameter
402856	Large Curl (LC)	80	110	7	Up to 0,035" guidewire
402857	Extra Large Curl (XLC)	80	110	7	Up to 0,035" guidewire



Large deflection for CS Cannulation

Small deflection for subselection of target branch vein

Separately Available Accessories

Model Number	Name	Pin Design	Working Length (cm)
402854	Bipolar Extension Cable	Shrouded 2 mm length pin	210
402855	Bipolar Extension Cable	Unshrouded 2 mm length pin	210

Customer Support: 46-8-474-4756

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use. Devices depicted may not be available in all countries. Check with your St. Jude Medical representative for product availability in your country. Unless otherwise noted, ™ indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2011 St. Jude Medical, Inc. All rights reserved.

CPS Venture™ Wire Control Catheter




Product Highlights

- Tip deflects up to 90° to steer any 0,014" guidewire through the most challenging venous anatomies
- Over-the-Wire design permits easy wire exchanges, if necessary

Ordering Information

Model Number	Working Length (cm)	Overall Length (cm)	Guidewire Compatibility (in)	Guide Catheter Compatibility (F)
1135-001	70	88	0,014	6



Customer Support: 46-8-474-4756

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use. Devices depicted may not be available in all countries. Check with your St. Jude Medical representative for product availability in your country. Unless otherwise noted, ™ indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2011 St. Jude Medical, Inc. All rights reserved.