

SJM Biocor™ Patch

Bovine Pericardium Repair Patch

Product Highlights

- Useful in a broad range of clinical indications
- Fixed with glutaraldehyde to increase stability and decrease antigenicity
- Available in multiple sizes to accommodate surgeon preference and need

Ordering Information

Contents: Repair Patch (1 unit per box)

Model/Reorder Number	Patch Size (cm)
B40	4 x 5
B40	5 x 7
B40	6 x 10
B40	10 x 12.5



Please review the Instructions for Use prior to using these devices for a complete listing of indications, contraindications, warnings, precautions, potential adverse events, and directions for use.

Product referenced is approved for CE Mark.

Unless otherwise noted, ™ indicates a registered trademark or service mark owned by, or licensed to, St. Jude Medical, Inc. or one of its subsidiaries. Biocor, ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL, LESS RISK, are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies.