

# RadiAnalyzer™ Xpress

## Instrumentation

### Product Highlights

- Used together with PressureWire™ Certus™ for real-time measurement of pressure (FFR), flow (CFR), microcirculatory resistance (IMR) and intravascular temperature
- Connects with any cath lab system
- Lightweight and compact monitor size allows bedside mounting for convenient access to key information
- High-speed connectivity to PC or printer

### Ordering Information

Contents: RadiAnalyzer Xpress (1 unit per box)

Note: 2 Monitor Cables and 1 Adapter Cable Required for Connection Need to Be Specified When Ordering.

Reorder Number	Description
12711	RadiAnalyzer Xpress



Product referenced is approved for CE Mark.