

Therapy™ Cool Path™ Ablation Catheter
 Irrigated 2 & 4 mm Tip
 Quadripolar
 7 F

Product Highlights

- Six equidistant ports on the distal tip
- Uni-directional steering
- Push/pull handle

Ordering Information

Contents: 7 F irrigated ablation catheter (1 unit per box)

Reorder Number	Description	Electrode Spacing (mm)	Tip Electrode (mm)	Band Electrode (mm)	Curve Type	Usable Length (cm)
83352	1304-CP-7-25-M-TC	2-5-2	4	2	Medium	110
83353	1304-CP-7-25-L-TC	2-5-2	4	2	Large	110
83354	1304-CP-7-25-XL-TC	2-5-2	4	2	X-Large	110
83358	1304-CP-7-25-M-TH	2-5-2	4	2	Medium	110
83359	1304-CP-7-25-L-TH	2-5-2	4	2	Large	110
83360	1304-CP-7-25-XL-TH	2-5-2	4	2	X-Large	110
83362	1304-CP-7-25-L-TE2-TC	2-5-2	2	2	Large	110
83376	1304-CP-7-25-FL-TC	2-5-2	4	2	Flutter	110
83382	1304-CP-7-25-L1-TC	2-5-2	4	2	Large1	110
83388	1304-CP-7-25-L1-TC	2	4	2	Large1	110


Required Catheter Connecting Cable – Page 165

SJM: 85641

Stockert: 85713, 85709

Medtronic: 85711, 85708

Biosense Webster: 85761



Therapy™ Cool Path™ Ablation Catheters

4 mm Tip Thermocouple

Quadripolar

7 F

Product Highlights

- All braided shaft
- Double mesh stainless steel braiding to the distal tip
- Six equidistant ports on the distal tip
- Uni-directional steering
- Push/pull handle

Ordering Information

Contents: 7 F irrigated ablation catheter (1 unit per box)


Reorder Number	Description	Electrode Spacing (mm)	Tip Electrode (mm)	Band Electrode (mm)	Curve Type	Usable Length (cm)
83502	1304-CP-7-25-L-AB	2-5-2	4	2	Large	110
83503	1304-CP-7-25-M-AB	2-5-2	4	2	Medium	110
83504	1304-CP-7-25-XL-AB	2-5-2	4	2	X-Large	110
83505	1304-CP-7-25-L1-AB	2-5-2	4	2	Large1	110
83506	1304-CP-7-25-FL-AB	2-5-2	4	2	Flutter	110

Required Catheter Connecting Cable – Page 165

SJM: 85641

Stockert: 85713, 85709

Medtronic: 85711, 85708



Therapy™ Cool Path™ Duo Ablation Catheters
 4 mm Tip Thermocouple
 Quadripolar
 7 F

Product Highlights

- All braided shaft
- Double mesh stainless steel braiding to the distal tip
- Two rows of six equidistant ports offset by 30 degrees
- Proximal and distal tip electrode irrigation
- Uni-directional steering
- Push/pull handle

Ordering Information

Contents: 7 F irrigated ablation catheter (1 unit per box)

Reorder Number	Description	Electrode Spacing (mm)	Tip Electrode (mm)	Band Electrode (mm)	Curve Type
83562	1304-CP2-7-1.5(5)2-M-TE4BE1AB	1.5-5-2	4	1	Medium
83563	1304-CP2-7-1.5(5)2-L-TE4BE1AB	1.5-5-2	4	1	Large
83564	1304-CP2-7-1.5(5)2-XL-TE4BE1AB	1.5-5-2	4	1	X-Large
83565	1304-CP2-7-1.5(5)2-L1-TE4BE1AB	1.5-5-2	4	1	Large1

Required Catheter Connecting Cable – Page 165

SJM: 85641, 85644

Stockert: 85713, 85709

Medtronic: 85711, 85708

