

Contact™ Therapy™ Cool Path™ Duo Ablation Catheter
 Irrigated 4 mm Tip Thermocouple
 Quadripolar

7 F

Product Highlights

- Compatible with EnSite™ Contact™ Technology
- Uni-directional steering with push/pull handle
- Proximal and distal tip electrode irrigation
- Two rows of six equidistant ports offset by 30 degrees
- Double mesh stainless steel braiding to the distal tip

Ordering Information

Contents: 7 F irrigated ablation catheter (1 unit per box)

Reorder Number	Description	Electrode Spacing (mm)	Tip Electrode (mm)	Band Electrode (mm)	Distal Reach	Usable Length (cm)
A088019	1304-CP2-7-1.5(5)2-M-TE4BE1AB-CN	1.5-5-2	4	1	Medium	110
A088020	1304-CP2-7-1.5(5)2-L-TE4BE1AB-CN	1.5-5-2	4	1	Large	110
A088022	1304-CP2-7-1.5(5)2-L1-TE4BE1AB-CN	1.5-5-2	4	1	L1	110

Required Catheter Connecting Cables – Page 167

Reorder Number

A402892

