

Safire™ Bi-directional Ablation Catheters

4 & 8 mm Tip Universal Temperature Monitoring

Quadripolar

7 F

Product Highlights

- Automatic steering lock reduces catheter manipulation
- Bi-directional tip deflection
- Thermocouple and Thermistor temperature monitoring capability
- ComfortGrip™ ergonomic handle

Ordering Information

Contents: 7 F Quadripolar catheter (1 unit per box)


Reorder Number	Electrode Spacing (mm)	Tip Electrode (mm)	Band Electrode (mm)	Distal Reach	Usable Length (cm)
402818	2-5-2	4	1	Small Sweep	115
402819	2-5-2	4	1	Medium Sweep	115
402820	2-5-2	4	1	Large Sweep	115
402821	2-5-2	4	1	Small Curl	115
402822	2-5-2	4	1	Medium Curl	115
402823	2-5-2	4	1	Large Curl	115
402824	2-5-2	8	1	Small Sweep	115
402825	2-5-2	8	1	Medium Sweep	115
402826	2-5-2	8	1	Large Sweep	115
402828	2-5-2	8	1	Medium Curl	115
402829	2-5-2	8	1	Large Curl	115

Required Catheter Extension Cables – Page 164

SJM: 402560, 402561

Stockert: 402566, 402567

EPT: 402558, 402559



Information contained herein for distribution outside of the U.S. only.

Electrophysiology International Catalogue

©2013 St. Jude Medical. All Rights Reserved.