

Therapy™ Ablation Catheters

4 mm Tip Thermocouple

5 F

Product Highlights

- Uni-directional steering
- Thermocouple temperature measurement
- Push/pull handle

Ordering Information

Contents: 7 F Quadripolar catheter (1 unit per box)


Reorder Number	Description	Electrode Spacing (mm)	Tip Electrode (mm)	Band Electrode (mm)	Curve Type
83351	1304-5-25-M-TE4BE1	2-5-2	4	1	Medium
83403	1304-5-2-S-TE4BE1	2	4	1	Small

Required Catheter Connecting Cable – Page 165

SJM: 85641

Stockert: 85713, 85709

Medtronic: 85711, 85708



Therapy™ Ablation Catheters

4 & 8 mm Tip Thermocouple

Quadripolar

7 F

Product Highlights

- Uni-directional steering
- Thermocouple temperature measurement
- Push/pull handle

Ordering Information

Contents: 7 F ablation catheter (1 unit per box)


Reorder Number	Description	Electrode Spacing (mm)	Tip Electrode (mm)	Band Electrode (mm)	Curve Type	Usable Length (cm)
83302	1304-7-25-M-TE8	2-5-2	8	2	Medium	110
83306	1304-7-25-L-TE8	2-5-2	8	2	Large	110
83309	1304-7-25-L-TE8A	2-5-2	8	2	Large	110
83311	1304-7-25-XL-TE8A	2-5-2	8	2	X-Large	110
83312	1304-7-25-XL-TE8	2-5-2	8	2	X-Large	110
83404	1304-7-2-M	2	4	2	Medium	110
83405	1304-7-25-M	2-5-2	4	2	Medium	110
83408	1304-7-25-L	2-5-2	4	2	Large	110
83411	1304-7-25-E	2-5-2	4	2	Extended Reach	110
83413	1304-7-2-F	2	4	2	Far Reach	110
83417	1304-7-25-XL	2-5-2	4	2	X-Large	110
83432	1304-7-25-S	2-5-2	4	2	Small	110

Required Catheter Connecting Cable – Page 165

SJM: 85641

Stockert: 85713, 85709

Medtronic: 85711, 85708



Therapy™ Ablation Catheters
 4 & 8 mm Tip Thermistor
 Quadripolar
 7 F

Product Highlights

- Uni-directional steering
- Thermistor temperature measurement
- Compatible with EPT-1000 XP generator
- Push/pull handle


Ordering Information

Contents: 7 F ablation catheter (1 unit per box)

Reorder Number	Description	Electrode Spacing (mm)	Tip Electrode (mm)	Band Electrode (mm)	Curve Type	Usable Length (cm)
83453	1304-7-25-S-TH	2-5-2	4	2	Small	110
83456	1304-7-25-M-TH	2-5-2	4	2	Medium	110
83459	1304-7-25-L-TH	2-5-2	4	2	Large	110
83462	1304-7-25-XL-TH	2-5-2	4	2	X-Large	110
83474	1304-7-25-M-TE8-TH	2-5-2	8	2	Medium	110
83477	1304-7-25-L-TE8-TH	2-5-2	8	2	Large	110
83481	1304-7-25-XL-TE8-TH	2-5-2	8	2	X-Large	110

Required Catheter Connecting Cable – Page 165

EPT: 85763



Therapy™ Dual-8™ Ablation Catheters

8 mm Tip Dual Thermocouple Quadripolar

7 F

Product Highlights

- Uni-directional steering
- 8 mm tip electrode with two thermocouples
- High power option up to 100 watts
- Push/pull handle

Ordering Information


Contents: 7 F ablation catheter (1 unit per box)

Reorder Number	Description	Electrode Spacing (mm)	Tip Electrode (mm)	Band Electrode (mm)	Curve Type	Usable Length (cm)
83421	1304-7-2-M-TE8TC2	2	8	2	Medium	110
83422	1304-7-25-M-TE8TC2	2-5-2	8	2	Medium	110
83424	1304-7-2-L-TE8TC2	2	8	2	Large	110
83425	1304-7-25-L-TE8TC2	2-5-2	8	2	Large	110
83426	1304-7-5-L-TE8TC2	5	8	2	Large	110
83427	1304-7-2-XL-TE8TC2	2	8	2	X-Large	110
83428	1304-7-25-XL-TE8TC2	2-5-2	8	2	X-Large	110
83307	1304-7-25-XL-TE8TC2A	2-5-2	8	2	X-Large	110
83308	1304-7-25-L-TE8TC2A	2-5-2	8	2	Large	110

Required Catheter Connecting Cable – Page 165

SJM: 85641

Stockert: 85739, 85643



Therapy™ Ablation Catheters
 4 & 8 mm Tip Thermocouple Triflex
 Quadripolar
 7 F

Product Highlights

- Uni-directional steering
- 4 mm tip with one thermocouple; 8 mm tip with two thermocouples
- Push/pull handle

Ordering Information

Contents: 7 F ablation catheter (1 unit per box)


Reorder Number	Description	Electrode Spacing (mm)	Tip Electrode (mm)	Band Electrode (mm)	Curve Type	Usable Length (cm)
83510	1304-7-25-M-TF	2-5-2	4	2	Medium	110
83513	1304-7-25-S-TF	2-5-2	4	2	Small	110
83516	1304-7-25-L-TF	2-5-2	4	2	Large	110
83542	1304-7-25-L-TE8TC2-TF	2-5-2	8	2	Large	110
83543	1304-7-25-XL-TE8TC2-TF	2-5-2	8	2	X-Large	110
83545	1304-7-25-M-TE8TC2-TF	2-5-2	8	2	Medium	110
83552	1304-7-25-M-TH-TF(90)	2-5-2	4	2	Medium	90

Required Catheter Connecting Cable – Page 165

SJM: 85641

Stockert: 85713, 85709

Medtronic: 85711, 85708



Therapy™ Ablation Catheters

4 & 8 mm Tip Thermocouple Triflex

Quadripolar

7 F

Product Highlights

- Bi-directional steering
- Thermocouple temperature measurement
- Push/pull handle

Ordering Information

Contents: 7 F ablation catheter (1 unit per box)

Reorder Number	Description	Electrode Spacing (mm)	Tip Electrode (mm)	Band Electrode (mm)	Curve Type
83701	1304-7-25-M-BD	2-5-2	4	2	Medium
83702	1304-7-25-L-BD	2-5-2	4	2	Large
83703	1304-7-25-S-BD	2-5-2	4	2	Small
83704	1304-7-25-L-TE8-BD	2-5-2	8	2	Large
83707	1304-7-25-M-TE8-BD	2-5-2	8	2	Medium
83708	1304-7-25-XL-TE8-BD	2-5-2	8	2	X-Large

Required Catheter Connecting Cable – Page 165

SJM: 85641

Stockert: 85713, 85709

Medtronic: 85711, 85708

