

# Livewire™ TC Ablation Catheters

## 4 mm Tip Universal Temperature Monitoring Quadripolar

### 7 F

### Product Highlights

- Bi-Directional tip deflection
- Thermocouple and Thermistor temperature monitoring capability
- Rotary dial designed for fine tip movements
- Extra lateral stability (XLS) incorporated at catheter tip
- Ergonomic handling

### Ordering Information

Contents: 7 F Quadripolar (XLS) catheter (1 unit per box)


Reorder Number	Electrode Spacing (mm)	Tip Electrode (mm)	Band Electrode (mm)	Distal Reach	Usable Length (cm)
402275	2-5-2	4	1	Med Sweep	115
402277	2-5-2	4	1	Med Curl	115
402278	2-5-2	4	1	Med Lrg Curl	115

#### Required Catheter Extension Cables – Page 163

SJM: 402543, 402544

Stockert: 402505

EPT: 402542



U.S. Patent Numbers 5,395,328 and 5,395,329

# Livewire™ TC Ablation Catheters

## 4 & 8 mm Tip Universal Temperature Monitoring

### Quadripolar

#### 7 F

### Product Highlights

- Thermocouple and Thermistor temperature monitoring capability
- Rotary dial designed for fine tip movements
- Extra lateral stability (XLS) incorporated at catheter tip
- Ergonomic handling
- Uni-directional tip deflection

### Ordering Information

Contents: 7 F Quadripolar (XLS) catheter (1 unit per box)


Reorder Number	Electrode Spacing (mm)	Tip Electrode (mm)	Band Electrode (mm)	Distal Reach	Usable Length (cm)
402112	2-5-2	4	1	Med Sweep	115
402113	2-5-2	4	1	Lrg Sweep	115
402114	2-5-2	4	1	Med Curl	115
402115	2-5-2	4	1	Lrg Curl	115
402193	2-5-2	8	1	Med Sweep	115
402194	2-5-2	8	1	Lrg Sweep	115
402195	2-5-2	8	1	Med Curl	115
402196	2-5-2	8	1	Lrg Curl	115

### Required Catheter Extension Cables – Page 163

SJM: 402543, 402544

Stockert: 402505

EPT: 402542



Information contained herein for distribution outside of the U.S. only.

Electrophysiology International Catalogue

©2013 St. Jude Medical. All Rights Reserved.

Livewire™ TC Uni-directional Ablation Catheters  
 4 & 8 mm Tip Thermocouple  
 Quadripolar  
 7 F

Product Highlights

- Thermocouple temperature monitoring capability
- Rotary dial designed for fine tip movements
- Ergonomic handling
- Uni-directional tip deflection

Ordering Information

Contents: 7 F Quadripolar (XLS) catheter (1 unit per box)

Reorder Number	Electrode Spacing (mm)	Tip Electrode (mm)	Distal Reach	Usable Length (cm)
402120	2-5-2	4	SmI Sweep	115
402121	2-5-2	4	Med Sweep	115
402122	2-5-2	4	Lrg Sweep	115
402124	2-5-2	4	Med Curl	115
402125	2-5-2	4	Lrg Curl	115
402104	2-5-2	8	Med Curl	115
402105	2-5-2	8	Lrg Curl	115

Required Catheter Extension Cables – Page 163

SJM: 402543, 402544

Stockert: 402505

