

## ViewFlex™ Xtra ICE Catheter

### Description

The ViewFlex Xtra ICE Catheter is a temporary intracardiac ultrasound catheter intended for use in patients to visualize cardiac structures and blood flow within the heart when connected to a ViewMate™ ultrasound console via the ViewFlex Catheter Interface Module.

### Product Highlights

- Single-handed, self-locking steering
- 9 F catheter size
- 10 F introducer
- 90 cm insertable length
- Four-way tip deflection, anterior/posterior, left/right
- 120 degree tip deflection
- 64 element phased array

### Ordering Information

Reorder Number	Description
100046963	ViewFlex Xtra ICE Catheter - Electrophysiology
100043720	ViewFlex PLUS ICE Catheter Interface Module for ViewMate Z Console
100038191	ViewFlex PLUS ICE Catheter Interface Module for ViewMate II System

