

# ViewMate™ Z Intracardiac Ultrasound Console

Powered By  
ZONARE Zone Sonography Technology



## EXCEPTIONAL, REAL-TIME VISUALIZATION FOR STREAMLINED PROCEDURES

The ViewMate Z Ultrasound Console is a full featured Intracardiac Echocardiography (ICE) system designed for exceptional, real-time image guidance and visualization of anatomical structures. The compact console delivers fast, high fidelity imaging performance to assist clinicians while performing complex electrophysiology and interventional procedures.

### Exceptional Visualization


- Real-time imaging streamlines procedures
- Premium image quality provides crisp, clear images of anatomical structures
- Software based ultrasound console allows for rapid upgrades providing the most up-to-date image quality

### Simplified Operation

- On-screen menus and displays help guide clinicians through the imaging process
- Presets and soft keys allow clinicians to save and retrieve settings effortlessly
- System recognition of transducer parameters optimize menu displays simplifying the learning curve
- Rapid boot-up time reduces time to begin procedures


### ZONARE Zone Sonography Technology

#### Conventional



*Line-by-line data acquisition limits image processing and formation.*

#### Zone Imaging



*Zone Sonography Technology decreases image acquisition time while improving image quality*



*Converts from a full featured cart based console into the Scan Engine, a 5½ lb. hand-held console without sacrificing performance.*



**ST. JUDE MEDICAL™**

MORE CONTROL. LESS RISK.

## Small footprint and lightweight system designed for effortless maneuverability and maximum versatility in tight or crowded spaces

- Single-handed release of 5½ lb. Scan Engine for portable scanning
- Rapid boot-up time reduces delay in starting procedures
- Multi-Transducer Port (MTP) - up to three transducers connected simultaneously
- Fully charged optional battery pack offers up to three hours of cordless operation
- 19" HD LCD with monitor mounted high-fidelity stereo speakers
- ViewFlex™ Catheter Interface Module
- Integrated front handle for quick vertical positioning and transport

### Ordering Information

Reorder Number	Product Description
100044000	ViewMate Z Console
100044271	ViewMate Z Console with battery pack
100042294	Transesophageal Transducer
100042295	Continuous Wave Transducer
100042741	Linear Transducer
100042293	Transthoracic Transducer
100043997	ViewMate Cart Battery Pack
100043998	Audio Video Extension Device
100043999	ECG Accessory Cable
100026714	ViewFlex™ PLUS ICE Catheter for ViewMate Z Console

ATRIAL FIBRILLATION

CARDIAC RHYTHM MANAGEMENT

CARDIOVASCULAR

NEUROMODULATION

#### Global Headquarters

One St. Jude Medical Drive  
St. Paul, Minnesota 55117  
USA  
+1 651 756 2000  
+1 651 756 3301 Fax

#### Atrial Fibrillation Division

One St. Jude Medical Drive  
St. Paul, Minnesota 55117  
USA  
+1 651 756 2000  
+1 651 756 3301 Fax

#### St. Jude Medical Coordination Center BVBA

The Corporate Village  
Da Vincilaan 11 Box F1  
1935 Zaventem, Belgium  
+32 2 774 68 11  
+32 2 772 83 84 Fax

#### St. Jude Medical Brasil, Ltda.

Rua Frei Caneca, 1380  
7º ao 9º andares  
01307-002 – São Paulo (SP)  
Brazil  
+55 11 5080 5400  
+55 11 5080 5423 Fax

#### St. Jude Medical (Hong Kong) Ltd.

Suite 1608, 16/F Exchange Tower  
33 Wang Chiu Road  
Kowloon Bay, Kowloon  
Hong Kong SAR  
+852 2996 7688  
+852 2956 0622 Fax

#### St. Jude Medical Japan Co., Ltd.

Shiodome City Center 15F  
1-5-2, Higashi-Shinbashi  
Minato-ku, Tokyo 105-7115  
Japan  
Tel: +81 3 6255 6370  
Fax: +81 3 6255 6371

[SJMprofessional.com](http://SJMprofessional.com)



**ST. JUDE MEDICAL™**  
MORE CONTROL. LESS RISK.

#### Rx Only

**Brief Summary:** Please review the Instructions for Use prior to using these devices for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

ZONARE Zone Sonography Technology is a registered or unregistered trademark of ZONARE Medical Systems, Inc.

ViewMate, ViewFlex, ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2011 St. Jude Medical. All Rights Reserved.

Item 100046964EN Rev B

