Myodex™

Bipolar, Steroid-eluting Epicardial Pacing Lead

Product Highlights

- Active-fixation, sutureless design
 - Full 3.5 mm helix penetration helps provide stable fixation
- Superb deliverability combined with exceptional stability and performance
- Low pacing thresholds with steroid elution
- Precise sensing and low polarisation
- Easy to implant with the FasTac[™] Lead Implant Tool



Ordering Information

Contents: Epicardial lead

B 4		I - I	
IVЛ	od		

Number	Insulation	Connector	Lengths (cm)
1084T	Silicone	IS-1 bipolar	25; 35; 54





Bipolar, Steroid-eluting Epicardial Lead

Product Specifications

DDELS	1084T (25 cm)	1084T (35 cm)	1084T (54 cm)
ectrode surface area	10 mm² cathode	10 mm ²	10 mm ²
	62 mm² anode	62 mm ²	62 mm ²
lix penetration depth	3,5 mm	3,5 mm	3,5 mm
mber of helix turns to implant	2,5	2,5	2,5
ad resistance	20 ohms cathode	27 ohms	41 ohms
	38 ohms anode	46 ohms	75 ohms
roducer length	27 cm	27 cm	27 cm
nneler length	27 cm	27 cm	27 cm
nnector type	IS-1 Bi	IS-1 Bi	IS-1 Bi
ectrode material	Titanium-nitride coated helix	Titanium-nitride coated	Titanium-nitride coated
	Platinum/Iridium	Platinum/Iridium	Platinum/Iridium
	Titanium-nitride coated anode plate	Titanium-nitride coated	Titanium-nitride coated
	titanium	titanium	titanium
nductor material	MP35N™ (multifilar coil)	MP35N™ (multifilar coil)	MP35N™ (multifilar coil)
sulation material	Silicone rubber	Silicone rubber	Silicone rubber
	(medical grade)	(medical grade)	(medical grade)
nnector pin material	316L stainless steel	316L stainless steel	316L stainless steel
eroid plug	<1 mg dexamethasone	<1 mg dexamethasone	<1 mg dexamethasone
	sodium phosphate	sodium phosphate	sodium phosphate
ture sleeve and pin cap material	Silicone rubber	Silicone rubber	Silicone rubber

ACCESSORIES PACKAGED WITH THE MYODEX LEAD

- 1 FasTac introducer
- 1 bidirectional tunneler tip
- 1 connector pin cap
- 1 slit suture sleeve (detached from lead)

