

St. Jude Medical Infinity™ IPG

Models 6660, 6661, 6662 and 6663



Compatibility

- St. Jude Medical™ Patient Controller Model 3875
- St. Jude Medical™ Clinician Programmer Model 3874

Hardware Specifications

St. Jude Medical Infinity™ IPG (Models 6660 and 6661*)

Specifications	Metric	Imperial
Height	5.55 cm	2.19 in
Length	4.95 cm	1.95 in
Thickness	1.34 cm	0.53 in
Weight	48.9 g	1.7 oz
Volume	30.4 cm ³	1.9 in ³
Battery	Carbon monofluoride/silver vanadium oxide cell	

St. Jude Medical Infinity™ IPG (Models 6662 and 6663*)

Specifications	Metric	Imperial
Height	6.68 cm	2.63 in
Length	5.02 cm	1.98 in
Thickness	1.35 cm	0.53 in
Weight	58.3 g	2.1 oz
Volume	38.6 cm ³	2.4 in ³
Battery	Carbon monofluoride/silver vanadium oxide cell	

Stimulation Performance

Amplitude	0-12.75 mA (0.05 to 1.00 mA steps)
Pulse Width Tonic	20-500 μs (10 μs steps)
Frequency	2-240 Hz (2 Hz steps)
Program Storage Capacity	16 programs

Accessories

Description	Model Number
Torque wrench	1101
Port plugs	1111
St. Jude Medical™ 8-channel adapter - B	2321 and 2326

Leads and Extensions Compatibility

6100-series	Traditional leads
6100-series	St. Jude Medical Infinity™ leads
6300-series**	Traditional extensions
6300-series**	St. Jude Medical Infinity™ extensions



*Models 6661 and 6663 are directly compatible with Medtronic extension series 37086

**Excludes model numbers 6315, 6316 and 6319

Rx Only

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

CE marking does not necessarily indicate regulatory approval status for all markets. Please refer to the instructions for use for a full listing of indications, contraindications, warnings and precautions.

Unless otherwise noted, TM indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL and the nine-squares symbol are trademarks and service marks of St. Jude Medical, Inc. and its related companies. © 2015 St. Jude Medical, Inc. All Rights Reserved.

St. Jude Medical Inc.

Global Headquarters
One St. Jude Medical Drive
St. Paul, MN 55117
USA
T +1 651 756 2000 | F +1 651 756 3301

St. Jude Medical S.C., Inc.

Americas Division
6300 Bee Cave Road
Bldg. Two, Suite 100
Austin, TX 78746
USA
T +1 512 286 4000 | F +1 512 732 2418

SJM Coordination Center BVBA

The Corporate Village
Da Vincilaan 11-Box F1
B-1935 Zaventem, Belgium
T +32 2 774 68 11 | F +32 2 772 83 84

St. Jude Medical Brasil Ltda.

Rua Itapeva, 538
5º ao 8º andares
01332-000 – São Paulo – SP
Brazil
T +55 11 5080 5400 | F +55 11 5080 5423

St. Jude Medical (Hong Kong) Ltd.

Suite 1608, 16/F Exchange Tower
33 Wang Chiu Road
Kowloon Bay, Kowloon
Hong Kong SAR
T +852 2996 7688 | F +852 2956 0622

St. Jude Medical Japan Co., Ltd.

Shiodome City Center 15F
1-5-2 Higashi Shinbashi, Minato-ku
Tokyo 105-7115
Japan
T +81 3 6255 6370 | F +81 3 6255 6371

St. Jude Medical Australia Pty, Ltd.

17 Orion Road
Lane Cove, NSW 2066
Australia
T +61 2 9936 1200 | F +61 2 9936 1222

SJMprofessional.com