



Needle Holders
with TC-Inserts



Driving your needles to perform




Regular	4-5
Fine	6-9
Fine Spring Style	10-11
Special	12-17
Short	18-19
With cutting function	20-21
Palm-held	22-23
Left-Hand	24-25
For Wire	26-33
Without TC-Inserts	34-37





TC-Insert Styles


Serration pitch



(S) **SAPPHIRE™**
Enhanced Needle Grip Surface




(S) **Smooth**
0 teeth PSI



(M) **Micro**
16000 teeth PSI
0.2 – 0.3 mm



(F) **Fine**
4000 teeth PSI
0.4 mm



(N) **Standard**
2600 teeth PSI
0.5 mm



(C) **Coarse**
1500 teeth PSI
0.65 mm



DeBakey



Mayo-Hegar



Micro-Vascular



Crile-Wood



DeBakey

Fig.	Art.-No.	length mm / "	Insert / Tip
A	24-0360*	180 / 7	ⓔ 1,5 mm
A	24-0361	200 / 8	ⓔ 1,5 mm
A	24-0362	230 / 9	ⓔ 1,5 mm
A	24-0363	260 / 10 ¼	ⓔ 1,5 mm
A	24-0364	300 / 11 ¾	ⓔ 1,5 mm

Titanium

Fig.	Art.-No.	length mm / "	Insert / Tip
A	24-0360.Ti	180 / 7	ⓔ 1,5 mm
A	24-0361.Ti	200 / 8	ⓔ 1,5 mm
A	24-0363.Ti	250 / 10	ⓔ 1,9 mm



A



C

Mayo-Hegar

Fig.	Art.-No.	length mm / "	Insert / Tip
C	24-0129	140 / 5 ½	Ⓝ 2,0 mm
C	24-0130	160 / 6 ¼	Ⓝ 2,5 mm
C	24-0131	180 / 7	Ⓝ 2,7 mm
C	24-0132*	200 / 8	Ⓝ 2,7 mm
C	24-0133	240 / 9 ½	Ⓝ 2,7 mm
C	24-0134	260 / 10 ¼	Ⓝ 2,7 mm
C	24-0135	300 / 11 ¾	Ⓝ 2,7 mm

Titanium

Fig.	Art.-No.	length mm / "	Insert / Tip
C	24-0130.Ti	160 / 6 ¼	Ⓝ 2,5 mm
C	24-0131.Ti	180 / 7	Ⓝ 2,7 mm
C	24-0132.Ti	200 / 8	Ⓝ 2,7 mm

Micro-Vascular

(DeBakey delicate)

Fig.	Art.-No.	length mm / "	Insert / Tip
B	24-0350*	150 / 6	Ⓜ 1,2 mm
B	24-0351	180 / 7	Ⓜ 1,2 mm
B	24-0352	200 / 8	Ⓜ 1,2 mm
B	24-0353	230 / 9	Ⓜ 1,2 mm



B




D

Crile-Wood


Fig.	Art.-No.	length mm / "	Insert / Tip
D	24-0310	150 / 6	ⓔ 2,0 mm
D	24-0311*	180 / 7	ⓔ 2,0 mm
D	24-0312	200 / 8	ⓔ 2,0 mm
D	24-0313	230 / 9	ⓔ 2,0 mm
D	24-0314	260 / 10 ¼	ⓔ 2,0 mm
D	24-0315	300 / 11 ¾	ⓔ 2,0 mm
D	24-0308	150 / 6	Ⓢ 1,7 mm

* bold = shown





see page 3



Sartorius



Bavaria



Geister-Vascular-Ryder



Geister-Vascular



Geister-Vascular Titanium



Geister-Vascular

Fig.	Art.-No.	length mm / "	Insert / Tip
C	24-0496	130 / 5	◆ 1,2 mm
C	24-0498	165 / 6 ½	◆ 1,2 mm
C	24-0500*	180 / 7	◆ 1,2 mm
C	24-0502	200 / 8	◆ 1,2 mm
C	24-0504	230 / 9	◆ 1,2 mm
C	24-0506	250 / 10	◆ 1,2 mm

Titanium

Fig.	Art.-No.	length mm / "	Insert / Tip
C	24-0530	125 / 5	◆ 1,2 mm
C	24-0532	145 / 5 ¾	◆ 1,2 mm
C	24-0534	175 / 7	◆ 1,2 mm
C	24-0536*	200 / 8	◆ 1,2 mm
C	24-0538	225 / 8 ¾	◆ 1,2 mm

For a secure but gentle grip of needles and sutures without damage in General Surgery, Cardiovascular Surgery, Plastic Surgery, Urology

Geister-Vascular-Ryder

Fig.	Art.-No.	length mm / "	Insert / Tip
D	24-0497	130 / 5	◆ 1,2 mm
D	24-0499	170 / 6 ¾	◆ 1,2 mm
D	24-0501*	180 / 7	◆ 1,2 mm
D	24-0503	190 / 7 ½	◆ 1,2 mm
D	24-0505	230 / 9	◆ 1,2 mm
D	24-0507	250 / 10	◆ 1,2 mm

Bavaria



A

Fig.	Art.-No.	length mm / "	Insert / Tip
A	24-0510*	230 / 9	Ⓝ 2,1 mm
A	24-0512	260 / 10 ¼	Ⓝ 2,1 mm

Ideally for larger needles and procedures demanding superior needle control.

Sarot



B

Fig.	Art.-No.	length mm / "	Insert / Tip
B	24-0331*	180 / 7	ⓔ 1,9 mm
B	24-0332	200 / 8	ⓔ 1,9 mm
B	24-0333	230 / 9	ⓔ 1,9 mm
B	24-0334	260 / 10 ¼	ⓔ 1,9 mm

* bold = shown



C



D





Mini-Ryder

Intracardiac

French-Eye (Ryder)

Eufrate-Pasqué

Neurosurgical

Lemmon



Neurosurgical

(French Eye (Ryder) delicate)



A

Fig.	Art.-No.	length mm / "	Insert / Tip
A	24-0346	130 / 5	ⓕ 1,7 mm
A	24-0348*	130 / 5	Ⓢ 1,7 mm
A	24-0347	150 / 6	ⓕ 1,7 mm
A	24-0349	150 / 6	Ⓢ 1,7 mm

French-Eye (Ryder)



B

Fig.	Art.-No.	length mm / "	Insert / Tip
B	24-0341	160 / 6 ¼	ⓕ 2,1 mm
B	24-0342	180 / 7	ⓕ 2,1 mm
B	24-0343*	200 / 8	ⓕ 2,1 mm
B	24-0344	240 / 9 ½	ⓕ 2,1 mm
B	24-0345	260 / 10 ¼	ⓕ 2,1 mm

Mini-Ryder

Fig.	Art.-No.	length mm / "	Insert / Tip
C	24-0410	130 / 5	Ⓜ 1,2 mm
C	24-0412	150 / 6	Ⓜ 1,2 mm
C	24-0413	180 / 7	Ⓜ 1,2 mm
C	24-0414*	200 / 8	Ⓜ 1,2 mm
C	24-0415	230 / 9	Ⓜ 1,2 mm

Titanium



C

Fig.	Art.-No.	length mm / "	Insert / Tip
C	24-0410.TI	130 / 5	◆ 1,2 mm
C	24-0412.TI	160 / 6 ¼	◆ 1,8 mm
C	24-0413.TI	180 / 7	◆ 1,8 mm
C	24-0414.TI	200 / 8	◆ 1,8 mm



D

Fig.	Art.-No.	length mm / "	Insert / Tip
D	24-0400	145 / 5 ¾	Ⓜ 1,6 mm
D	24-0401*	170 / 6 ¾	Ⓜ 1,6 mm
D	24-0402	195 / 7 ¾	Ⓜ 1,6 mm

Eufrate-Pasqué



E

Fig.	Art.-No.	length mm / "	Insert / Tip
E	24-0444	140 / 5 ½	Ⓜ 1,2 mm
E	24-0445	150 / 6	Ⓜ 1,2 mm
E	24-0446*	180 / 7	Ⓜ 1,2 mm
E	24-0447	200 / 8	Ⓜ 1,2 mm
E	24-0439	140 / 5 ½	Ⓢ 1,2 mm
E	24-0440	150 / 6	Ⓢ 1,2 mm
E	24-0441	180 / 7	Ⓢ 1,2 mm
E	24-0442	200 / 8	Ⓢ 1,2 mm

Intracardiac



F

Fig.	Art.-No.	length mm / "	Insert / Tip
F	24-0322*	180 / 7	ⓕ 1,7 mm
F	24-0323	200 / 8	ⓕ 1,7 mm
F	24-0324	230 / 9	ⓕ 1,7 mm
F	24-0325	260 / 10 ¼	ⓕ 1,7 mm

* bold = shown



N



F



M



S



see page 3

Giannini


Giannini Heavy

Castroviejo

Diethrich Ryder

Castroviejo / Hepp-Scheidel

Fairgrip™ultra





Giannini Heavy

with ratchet

Fig.	Art.-No.	length mm / "	Insert / Tip
A	24-0258*	260 / 10 1/4	Ⓜ 1,9 mm



Giannini

with ratchet

Fig.	Art.-No.	length mm / "	Insert / Tip
B	24-0255*	240 / 9 1/2	Ⓜ 1,3 mm



Castroviejo / Hepp-Scheidel

with ratchet

Fig.	Art.-No.	length mm / "	Insert / Tip
C	24-0250	130 / 5	Ⓢ 1,2 mm
C	24-0252	130 / 5	Ⓜ 1,2 mm
C	03-0270*	180 / 7	Ⓜ 1,2 mm



Diethrich Ryder

with ratchet

Fig.	Art.-No.	length mm / "	Insert / Tip
D	24-0420	150 / 6	Ⓜ 1,2 mm
D	24-0421	180 / 7	Ⓜ 1,2 mm
D	24-0422	200 / 8	Ⓜ 1,2 mm
D	24-0423*	230 / 9	Ⓜ 1,2 mm
D	24-0424	250 / 10	Ⓜ 1,4 mm



Fairgrip™ ultra

with ratchet

Fig.	Art.-No.	length mm / "	Insert / Tip
E	24-1140	180 / 7	Ⓜ 1,2 mm
E	24-1150	210 / 8 1/4	Ⓜ 1,2 mm
E	24-1160	230 / 9	Ⓜ 1,2 mm
F	24-1141*	180 / 7	Ⓜ 1,2 mm
F	24-1151	210 / 8 1/4	Ⓜ 1,2 mm
F	24-1161	230 / 9	Ⓜ 1,2 mm




Castroviejo

with ratchet

Fig.	Art.-No.	length mm / "	Insert / Tip
G	24-0200	145 / 5 3/4	Ⓢ 1,2 mm
G	24-0220	145 / 5 3/4	Ⓜ 1,2 mm
H	24-0201	145 / 5 3/4	Ⓢ 1,2 mm
H	24-0221	145 / 5 3/4	Ⓜ 1,2 mm
G	24-0204	180 / 7	Ⓢ 1,2 mm
G	24-0224	180 / 7	Ⓜ 1,2 mm
H	24-0205	180 / 7	Ⓢ 1,2 mm
H	24-0225	180 / 7	Ⓜ 1,2 mm
G	24-0208	220 / 8 3/4	Ⓢ 1,2 mm
G	24-0228*	220 / 8 3/4	Ⓜ 1,2 mm
H	24-0209	220 / 8 3/4	Ⓢ 1,2 mm
H	24-0229	220 / 8 3/4	Ⓜ 1,2 mm

* bold = shown





Masson




Wangensteen



Mod. DeBakey



Hegar



Crile-Wood





Masson

Fig.	Art.-No.	length mm / "	Insert / Tip
C	24-0830*	265 / 10 1/2	Ⓝ 3,2 mm



Hegar

Fig.	Art.-No.	length mm / "	Insert / Tip
A	24-0150*	200 / 8	Ⓝ 3,5 mm



Wangensteen

Fig.	Art.-No.	length mm / "	Insert / Tip
D	24-0379*	270 / 10 1/2	Ⓝ 2,4 mm



Crile-Wood

with extended ring

Fig.	Art.-No.	length mm / "	Insert / Tip
B	24-0762	200 / 8	Ⓞ 2,0 mm
B	24-0764*	220 / 8 3/4	Ⓞ 2,0 mm
B	24-0766	240 / 9 1/2	Ⓞ 2,0 mm



Mod. DeBakey

with extended ring

Fig.	Art.-No.	length mm / "	Insert / Tip
E	24-0761	170 / 6 3/4	Ⓞ 1,7 mm
E	24-0763*	200 / 8	Ⓞ 1,7 mm
E	24-0767	250 / 10	Ⓞ 1,7 mm

* bold = shown







Bozemann-Finochietto



Crile-Wood



Heaney



Finochietto



Heaney

Fig.	Art.-No.	length mm / "	Insert / Tip
A	24-0145	210 / 8 ¼	ⓕ 2,6 mm
A	24-0146*	240 / 9 ½	ⓕ 2,6 mm
A	24-0147	270 / 10 ½	ⓕ 2,6 mm



Finochietto

Fig.	Art.-No.	length mm / "	Insert / Tip
B	24-0148*	210 / 8 ¼	ⓕ 2,5 mm
B	24-0149	270 / 10 ½	ⓕ 2,5 mm



Bozemmann-Finochietto

Fig.	Art.-No.	length mm / "	Insert / Tip
C	24-0158*	200 / 8	ⓕ 2,5 mm
C	24-0159	270 / 10 ½	ⓕ 2,5 mm



Crile-Wood

curved

Fig.	Art.-No.	length mm / "	Insert / Tip
D	24-0306	130 / 5	ⓕ 1,8 mm
D	24-0316	150 / 6	ⓕ 1,8 mm
D	24-0317	180 / 7	ⓕ 1,8 mm
D	24-0318*	200 / 8	ⓕ 1,8 mm
D	24-0319	230 / 9	ⓕ 2,0 mm

* bold = shown






Bozemann (Wertheim)

Masson-Lüthy

Stratte

Stratte reverse





Stratte reverse

Fig.	Art.-No.	length mm / "	Insert / Tip
C	24-0373	205 / 8	Ⓝ 2,7 mm
C	24-0371*	235 / 9 ¼	Ⓝ 2,7 mm




Masson-Lüthy

Fig.	Art.-No.	length mm / "	Insert / Tip
A	24-0156*	200 / 8	Ⓝ 2,5 mm
A	24-0157	260 / 10 ¼	Ⓝ 2,5 mm



Bozemann (Wertheim)

Fig.	Art.-No.	length mm / "	Insert / Tip
B	24-0140*	200 / 8	Ⓝ 2,5 mm
B	24-0141	230 / 9	Ⓝ 2,5 mm
B	24-0142	260 / 10 ¼	Ⓝ 2,5 mm
B	24-0143	290 / 11 ½	Ⓝ 2,5 mm



Stratte

Fig.	Art.-No.	length mm / "	Insert / Tip
D	24-0372	205 / 8	Ⓝ 2,7 mm
D	24-0370*	235 / 9 ¼	Ⓝ 2,7 mm

* bold = shown





see page 3

Baumgartner




Baby Crile-Wood



Halsey



Webster



Swedish



Derf





Baumgartner

Fig.	Art.-No.	length mm / "	Insert / Tip
A	24-0302*	130 / 5	Ⓝ 2,6 mm



Derf

Fig.	Art.-No.	length mm / "	Insert / Tip
B	24-0305*	120 / 4 ¾	ⓔ 1,4 mm



Swedish

Fig.	Art.-No.	length mm / "	Insert / Tip
C	24-0490*	140 / 5 ½	ⓔ 1,7 mm

* bold = shown



Webster

Fig.	Art.-No.	length mm / "	Insert / Tip
D	24-0292*	120 / 4 ¾	Ⓢ 1,4 mm
D	24-0291	120 / 4 ¾	Ⓜ 1,4 mm



Halsey

Fig.	Art.-No.	length mm / "	Insert / Tip
E	24-0301	130 / 5	ⓕ 1,5 mm
Titanium			
E	24-0301.Ti	130 / 5	ⓕ 1,8 mm



Baby Crile-Wood

Fig.	Art.-No.	length mm / "	Insert / Tip
F	24-0307	150 / 6	ⓕ 1,7 mm
Titanium			
F	24-0307.Ti	160 / 6 ¼	ⓕ 1,7 mm





Universal TC Wire Cutter

Kilner

Olsen-Hegar mini

Gillies

Neivert (Offset Ring)



Universal TC Wire Cutter

Fig. A Art.-No. **26-0140*** length mm / " 120 / 4 ¾

A



Olsen-Hegar mini with cutting edge

Fig. B Art.-No. **24-0160*** length mm / " 120 / 4 ¾ Insert / Tip (F) 1,8 mm

B



Neivert (Offset Ring)

Fig. C Art.-No. **24-0304*** length mm / " 130 / 5 Insert / Tip (M) 1,7 mm

C



Kilner

Fig. D Art.-No. **24-0303*** length mm / " 140 / 5 ½ Insert / Tip (E) 1,6 mm

D




Gillies with cutting edge

Fig. E Art.-No. **24-0167*** length mm / " 160 / 6 ¼ Insert / Tip (F) 2,2 mm

E

* bold = shown





Mathieu

Mathieu
Internal ratchet

Zweifel



Zweifel

Fig.	Art.-No.	length mm / "	Insert / Tip
A	24-0050*	200 / 8	Ⓝ 2,9 mm



Mathieu

Fig.	Art.-No.	length mm / "	Insert / Tip
B	24-0040	140 / 5 1/2	ⓕ 2,4 mm
B	24-0041	170 / 6 3/4	Ⓝ 2,4 mm
B	24-0042*	200 / 8	Ⓝ 2,4 mm



Mathieu

Internal ratchet


Fig.	Art.-No.	length mm / "	Insert / Tip
C	24-0046	170 / 6 3/4	Ⓝ 2,9 mm
C	24-0047*	200 / 8	Ⓝ 2,9 mm
C	24-0048	230 / 9	ⓕ 2,9 mm

* bold = shown





Left-hand French-Eye (Ryder)



Left-hand Mayo-Hegar



Left-hand DeBakey



Left-hand Baumgartner





Left-hand Mayo-Hegar

Fig.	Art.-No.	length mm / "	Insert / Tip
A	24-7530.LH	160 / 6 ¼	ⓕ 2,3 mm
A	24-7531.LH	190 / 7 ½	Ⓝ 2,3 mm
A	24-7532.LH*	200 / 8	Ⓝ 2,5 mm

A

A



Left-hand Baumgartner

Fig.	Art.-No.	length mm / "	Insert / Tip
C	24-7582.LH*	130 / 5	ⓕ 2,2 mm

C

Fig.

Art.-No.

length mm / "

Insert / Tip

C

24-7582.LH*

130 / 5

ⓕ 2,2 mm



Left-hand French-Eye (Ryder)

Fig.	Art.-No.	length mm / "	Insert / Tip
B	24-7560.LH	160 / 6 ¼	ⓕ 1,9 mm
B	24-7561.LH*	180 / 7	ⓕ 1,9 mm

B

B



Left-hand DeBakey

Fig.	Art.-No.	length mm / "	Insert / Tip
D	24-7550.LH	180 / 7	ⓕ 2,0 mm
D	24-7552.LH*	230 / 9	ⓕ 2,0 mm
D	24-7553.LH	260 / 10 ¼	ⓕ 2,0 mm
D	24-7554.LH	300 / 11 ¾	ⓕ 2,0 mm

D

Fig.

Art.-No.

length mm / "

Insert / Tip

D

24-7550.LH

180 / 7

ⓕ 2,0 mm

D

24-7552.LH*

230 / 9

ⓕ 2,0 mm

D

24-7553.LH

260 / 10 ¼

ⓕ 2,0 mm

D

24-7554.LH

300 / 11 ¾

ⓕ 2,0 mm

* bold = shown





see page 3



Berry



Rubio



Heavy Berry



A

Berry

Fig.	Art.-No.	length mm / "	Insert / Tip
A	24-0119*	180 / 7	Ⓝ 5,3 mm



B

Heavy Berry

Wire Twister with TC-Inserts

Fig.	Art.-No.	length mm / "	Insert / Tip
B	24-0127*	200 / 8	Ⓧ 6,0 mm




C

Rubio

Mini Wire Twister with TC-Inserts

Fig.	Art.-No.	length mm / "	Insert / Tip
C	24-0128*	130 / 5	Ⓧ 4,8 mm

* bold = shown





Wire Twister



Mathieu-Style




Automatic




Mathieu-Style

Fig.	Art.-No.	length mm / "	Insert / Tip
A	24-0980	170 / 6 3/4	Ⓝ 4,0 mm
B	24-0981*	200 / 8	Ⓝ 4,0 mm



Wire Twister


Fig.	Art.-No.	length mm / "	Insert / Tip
C	24-0130	160 / 6 1/4	Ⓝ 2,5 mm
D	24-0125*	155 / 6	Ⓝ 4,0 mm
D	24-0126	180 / 7	Ⓝ 4,0 mm
E	24-0120	155 / 6	Ⓝ 3,0 mm
E	24-0121	180 / 7	Ⓝ 3,0 mm



Automatic


Fig.	Art.-No.	length mm / "	Insert / Tip
F	24-0110*	160 / 6 1/4	Ⓝ 2,5 mm
G	24-0105	155 / 6	Ⓝ 4,0 mm
H	24-0100	155 / 6	Ⓝ 3,0 mm
I	24-0115	180 / 7	Ⓝ 5,3 mm

* bold = shown





see page 3



Wire Cutting Pliers

Wire Pliers

Wire-Twister

Wire Holding Forceps


Izzat Wire-Grip



A

Wire Pliers

Fig.	Art.-No.	length mm / "
A	14-0670*	170 / 6 ¾



B

C

D

Wire-Twister

Fig.	Art.-No.	length mm / "	Insert / Tip
B	24-0956	180 / 7	Ⓢ 2,0 mm
C	24-0955	180 / 7	Ⓢ 4,0 mm
D	24-0950*	180 / 7	Ⓢ 6,0 mm

* bold = shown

Wire Holding Forceps

Art.-No.	length mm / "
14-0700*	130 / 5


Wire Cutting Pliers

"removeable handles"

Art.-No.	length mm / "	Wire max. Ø
24-0930*	570 / 22 ½	6,35 mm

Izzat Wire-Grip

Art.-No.	Wire max. Ø
14-0710*	1,3 mm





Maxi Cutter




Wire Bull



Midi Cutter



Power Cutter



Wire Plier




A

1/2

Wire Bull

Fig.	Art.-No.	length mm / "	Wire max. Ø
A	24-0942	250 / 10	3,0 mm

Plastic inserts hold wire piece after being out



B

1/4

Wire Pliers


Fig.	Art.-No.	length mm / "
B	14-0675*	160 / 6 1/4

Midi Cutter

Art.-No.	length mm / "	Wire max. Ø
24-0906	180 / 7	1,5 mm
24-0905*	200 / 8	1,5 mm

Maxi Cutter

Art.-No.	length mm / "	Wire max. Ø
24-0910	175 / 7	1,7 mm
24-0920*	220 / 8 3/4	2,2 mm



C

1/2

Power Cutter

also recommended for plates

Fig.	Art.-No.	length mm / "	Wire max. Ø
C	24-0940*	230 / 9	3,5 mm

* bold = shown



Crile-Wood



Mayo-Hegar

Adson



Crile-Wood

Fig.	Art.-No.	length mm / "	Tip
A	24-8131*	150 / 6	2,0 mm



Mayo-Hegar

Fig.	Art.-No.	length mm / "	Tip
B	24-8142	160 / 6 ¼	2,7 mm
B	24-8143*	190 / 7 ½	2,7 mm
B	24-8144	200 / 8	3,0 mm
B	24-8145	240 / 9 ½	3,0 mm
B	24-8146	270 / 10 ½	3,0 mm



Adson

Fig.	Art.-No.	length mm / "	Tip
C	24-8160*	180 / 7	3,0 mm

* bold = shown



Webster



Mathieu



Baumgartner



Derf



A

Derf

Fig.	Art.-No.	length mm / "	Tip
A	24-8101*	120 / 4 ¾	2,0 mm



B

Webster

plain

Fig.	Art.-No.	length mm / "	Tip
B	24-8111*	130 / 5	2,2 mm



C

Baumgartner

Fig.	Art.-No.	length mm / "	Tip
C	24-8115*	130 / 5	2,5 mm



D



Mathieu

Fig.	Art.-No.	length mm / "	Tip
D	24-8170	140 / 5 ½	2,4 mm
D	24-8171*	170 / 6 ¾	2,4 mm
D	24-8172	200 / 8	2,4 mm

* bold = shown



We are focused on the fields „Cardio“ and „Neuro“ where we want to offer our valued clients a complete range of products and services. Our products are clearly described and illustrated in well-structured sales documents. A small selection of various catalogs we show you below.







THE NEW SET-CONFIGURATOR

In our new Set-Configurator you can find over 9000 special surgical instruments. Search systematically for specific product codes, keywords or reference numbers or click through the various index areas. Using the bookmark function, you can generate your individual product list, which you can print or save or fast and easily request a quote.

– www.geister.com –





GEISTER Medizintechnik GmbH

Föhrenstrasse 2 | D-78532 Tuttlingen | Germany

☎ +49 7461/966240 | 📠 +49 7461/9662422

www.geister.com | info@geister.com