









GEISTER®

SURGIMOTOR™

Our innovative new SURGIMOTOR[™] primary sternum saws with the latest Li-Ion battery technology feature a powerful engine integrated along with their electronics into the battery packs, which will be removed during sterilization.

Therefore an efficient cleaning and lubrication of the machines and their mechanics become possible and the lifetime can be increased. The compact size and low weight provide for comfortable handling. The sternum guard can be turned in 4 positions with 90° increments without tools.

SURGIMOTOR™ Sternum Saws with Li-Ion Power Technology

Art.	No.	Description	
14-6	5700	Sternum Saw reciprocating	
		Dimensions (without attachment)	217 x 93 x 197 mm
		Weight with PowerPack	1750 g
		Continuous adjustable rotation speed	0 – 6000 rpm (Step I), 0 – 6500 rpm (Step II)
14-6	5710	Sternum Saw Redo oscillating	
		Dimensions (without attachment)	163 x 93 x 202 mm
		Weight with PowerPack	1760 g
		Continuous adjustable rotation speed	0 – 9000 rpm (Step I), 0 – 11000 rpm (Step II)

Saw Blades for 14-6700

non-sterile

Art.-No.Description14-6705NSternum Saw Blade, 32 x 1.05 mm (0.8)

Saw Blades for 14-6710

non-sterile

ArtNo.	Description
14-6714N	Sternum Saw Blade oscillating, 35 x 20 mm
14-6715N	Sternum Saw Blade oscillating, 50 x 40 mm
14-6716N	Sternum Saw Blade oscillating, 35 x 75 mm





SURGIMOTOR™ Sternum Saws



SURGIMOTOR™ Tray

ArtNo.	Description
14-6740	Tray with Lid for 2 Sternum Saws and Sterile Shield
14-6740.BO	Tray without Lid
14-6740.LD	Lid for Tray
14-6741.BO	Tray without Lid for 1 Sterum Saw and Sterile Shield

SURGIMOTOR[™] Power Sources

ArtNo.	Description
14-6720	Li-Ion Power Pack
14-6725	Power Pack Charging Bay
14-6726	4 Power Pack Charging Bay

SURGIMOTOR™ Power Cable

Description
Cable for Charger Australia-plug
Cable for Charger EU-plug
Cable for Charger UK-plug
Cable for Charger USA-plug
POAG-Cable for Charger



Accessories for SURGIMOTOR™





For Trays and Containers see our special catalog.

Accessories for SURGIMOTOR™

ArtNo.	Description
14-6702	Sternum Saw Guard
14-6722	Sterile Shield
14-6728	Lubrication Stand
14-6729	Lubrication Oil Spray
14-6730	Spray Adapter, distal end, yellow
14-6731	Spray Adapter, center, blue
14-6732	Lubrication Adapter, universal, lilac
14-6735	Cleaning Brushes Set
14-6736	Flushing Adapter, Sternum Saw, red
14-6737	Flushing Adapter, Sternum Saw, green
14-6738	3mach Cleaning Bridge
14-6739.BL	3mach Cleaning Bridge, Connector Belimed
14-6739.GE	3mach Cleaning Bridge,
	Connector Maquet/MediKomp/Getinge
14-6739.ML	3mach Cleaning Bridge,
	Connector Miele/Stielco/Webeco

GEISTER[®]



GEISTER Medizintechnik GmbH

Föhrenstrasse 2 | D-78532 Tuttlingen | Germany ☎ +49 7461/966240 | ♣ +49 7461/9662422 www.geister.com | info@geister.com