AAJUOSAVOIGHAO INJMJDANAM MHIYHY DAIQHAD NOTIALJIYAH JAIYIA

xe1 2200 0002 208+

Kowloon Bay, Kowloon

Suite 1608, 16/F Exchange Tower

St. Jude Medical (Hong Kong) Ltd.

33 Wang Chiu Road

+1 651 756 4466 Fax

St. Paul, Minnesota 55117

0244 327 126 4470

+852 2996 7688

AA2 gnoX gnoH

ΝΟΠΑΟΜΟΡΟΓΑΠΟΝ

XE1 9866 2046 5 18+

3-1-30, Minami-Aoyana

St. Jude Medical Japan Co., Ltd.

+81 🖁 3453 6420

105 JOL 0085

ueder

Minato-ku

377 East County Road B Cardiovascular Division

ASU

X64 1055 861 168 1+ +1 651 756 2000 ASU St. Paul, Minnesota 55117 One St. Jude Medical Drive

ABVB retrestion Center BVBA St. Jude Medical

Global Headquarters

+32 2 / / 2 83 84 Fax +32 2 774 68 11 muiglað metnevsZ 28ventem First and I The Corporate Village

moo.lenoizesfonal.com

Please review the Instructions for Use prior to using these devices for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use. Rx Only

Products referenced are approved for CE Mark.

and its related companies. @2011 St. Jude Medical. All rights reserved. Engage, RadiStop, ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK: are registered and unregistered trademarks and service marks of St. Jude Medical, Inc.

X61 5242 0800 11 66+

Rua Frei Caneca, 1380

+1 512 732 2418 Fax +1 512 732 7400

Austin, Texas 78746

001 stiu2 ,owT gnibliu8

6300 Bee Cave Road

noisivi**D**.2.U

01307-002 - São Paulo (SP)

St. Jude Medical Brasil Ltda.

+22 JJ 2080 2400

7° ao 9° andares


Brazil

ASU

Item 100052243 IPN 1384-11







ACCESS WITH CLOSURE IN MIND.

Engage[™] introducers from St. Jude Medical set the stage for successful patient outcomes, with features designed to provide reliable vascular access and reduce trauma, which have the potential to result in better closure. Our focus on the entire interventional procedure, combined with state-of-the-art manufacturing technology and regoraustesting, help ensure we deliver on our promise to provide only the best possible product. For efficient, consistent hemostasis introducer performance, choose Engage.

Thoughtfully designed, rigorously tested for efficient and consistent performance

Flexible suture ring

Physician-driven features improve control, reduce trauma

- The unique rotating strain relief and attached flexible suture ring provide sheath kink resistance and convenient suture placement
- A shorter dilator and lengthened tapers provide more control, helping reduce the insertion force required for access
- Smooth tapers and tip transitions allow consistent, less traumatic vascular access across a range of patient and procedure types

Improved access may lead to optimal closure

- Precise tipping enhances procedural control, which has the potential to result in better closure
- Smooth sheath and dilator transitions help reduce vascular trauma
- The combination of transition features helps enhance the entire interventional procedure





Design innovations enhance efficiency

- Seal balances pass-through friction, maintaining reliable hemostasis using devices up to 9 F. The seal is maintained even if no device is inserted
- Streamlined hub is designed to decrease internal clot formation

St. Jude Medical is committed to putting more control into the hands of physicians, thereby reducing risk. As part of a comprehensive portfolio that delivers on this promise, the Engage family of femoral and radial introducers sets the stage for successful interventional procedures by offering reliable vascular access across a broad range of patients and procedures, while enhancing efficiency and helping minimize procedural risk.

Ordering Information

Engage Introducer (Femoral)

Reorder Number	French Size (ACT)	Sheath Length	Guidewire Diameter (in / mm)
C408500	4 F	12 cm	0.035 / 0.889
C408501	5 F	12 cm	0.035 / 0.889
C408502	6 F	12 cm	0.035 / 0.889
C408503	7 F	12 cm	0.035 / 0.889
C408504	8 F	12 cm	0.035 / 0.889
C408505	9 F	12 cm	0.035 / 0.889
C408511	4 F	25 cm	0.035 / 0.889
C408512	5 F	25 cm	0.035 / 0.889
C408513	6 F	25 cm	0.035 / 0.889
C408514	7 F	25 cm	0.035 / 0.889
C408515	8 F	25 cm	0.035 / 0.889
C408516	9 F	25 cm	0.035 / 0.889

	Reorder Number	French Size (ACT)	Sheath Length	Guidewire Diameter (in / mm)
	C408506	5 F	12 cm	0.038/0.965
	C408507	6 F	12 cm	0.038 (0.965
	C408508	7 F	12 cm	0.038/0.965
	C408509	8 F	12 cm	0.038 / 0.965
	C408510	9 F	12 cm	0.038 / 0.965
	C408517	51	25 cm	0.038 / 0.965
	C408518	6 F	25 cm	0.038 / 0.965
	C408519	7.7	25 cm	0.038 / 0.965
	0408520	8 F	25 cm	0.038 / 0.965
6	C 408521	9 F	25 cm	0.038/0.965

Accessories (Femoral)

Kit Contents (10 units per box)

- Introducer sheath with hemostasis valve and side port
- Dilator
- 50 cm double distal guidewire with "J" and straight ends

- 80 cm double distal guidewire for 25 cm sheath lengths

Engage TR Introducer (Radial)

Reorder Number	French Size (ACT)	Wire/Needle Configuration	Sheath Length	Guidewire Diameter (in / mm)
C408541	4 F	SS/standard	8.5 cm	0.025 / 0.635
C408543	5 F	SS/standard	8.5 cm	0.025 / 0.635
C408545	6 F	SS/standard	8.5 cm	0.025 / 0.635
C408525	4 F	SS/standard	12 cm	0.025 / 0.635
C408526	5 F	SS/standard	12 cm	0.025 / 0.635
C408527	6 F	SS/standard	12 cm	0.025 / 0.635
C408528	4 F	SS/standard	25 cm	0.025 / 0.635
C408529	5 F	SS/standard	25 cm	0.025 / 0.635
C408530	6 F	SS/standard	25 cm	0.025 / 0.635

Reorder Number	French Size (ACT)	Wire/Nsedle Configuration	Sheath Length	Guidewire Diameter (in / mm)
C408540	4 F	Conted/cannula	8.5 cm	0.025 / 0.635
C408542	5 F	Coated/cannula	8.5 cm	0.025 / 0.635
C408544	6 F	Coated/cannula	8.5 cm	0.025 / 0.635
C408534	4 F	Coated/cannula	12 cm	0.025 / 0.635
C408535	5 F	Coated/cannula	12 cm	0.025 / 0.635
C408536	SF	Coated/cannula	12 cm	0.025 / 0.635
C-108537	4 F	Coated/cannula	25 cm	0.025 / 0.635
C408538	5 F	Coated/cannula	25 cm	0.025 / 0.635
C+08539	6 F	Coated/cannula	25 cm	0.025 / 0.635

Accessories (Radial)

Kit Contents (5 units per box)

- Introducer sheath with hemostasis valve and side port

- Dilator

- 50 cm (80 cm for 25 cm sheath lengths) single distal guidewire (stainle's steel or PTF5-coated) with floppy straight end 20 gauge, 3.8 cm standard stainless steel needle or 20 gauge 3.2 cm cannula needle