

GEISTER

Clean
solutions

Mod. Ankeney-type
Sternal Retractor

GEISTER

The better way to




Modified Ankeny-type Sternal Retractor



Clean solutions:

One of the most popular large sternal retractors for full sternotomy has been further advanced with quick-removable and exchangeable suture organizing springs to meet the requirements of modern CSSD units.



Modified Ankeny-type Sternal Retractor self-retaining with rack and pinion action. Complete with six detachable swivel blades, and one pair of quick-removable springs on each arm. One pack of replacement springs 29-1300.FE included.

29-1300

**Mod. Ankeny-type Sternal retractor with swivel blades, standard (6)
Blade Depth 32 mm
Blade Width 35 mm
Total Opening 178 mm**

29-1300.42

**Extra swivel blades, deep (6)
Blade Depth 42 mm
Blade Width 35 mm**

29-1300.52

**Extra swivel blades, extra deep (6)
Blade Depth 52 mm
Blade Width 35 mm**

29-1300.FE

**Replacement springs (pack of 4)
(not sterile)**



The better way to operate™

Geister Medizintechnik GmbH
Föhrenstrasse 2
D-78532 Tuttlingen/Germany
Phone +49 (0) 7461 – 966 240
Telefax +49 (0) 7461 – 966 2422
E-mail info@geister.com
Internet www.geister.com

Distributor