# Integra®

Surgical Illumination and Visualization Systems

Limit uncertainty with maximum light to visualize critical anatomy





# Table of Contents

LED Headlight System	3
Parts and Accessories	5
MLX 300	
UltraLite® Pro Headlight	7
UltraLite <sup>®</sup> Plus Headlight	3
DLX UltraLite <sup>®</sup> Pro Camera	
Parts and Accessories10	

# LED Headlight System

Limit uncertainty with the freedom of cool, bright light. Integra introduces the superior surgical LED headlight that delivers the cool brightness and color you need with increased freedom and lower cost of ownership.



### Cool, bright LED light

- Homogenous illumination across a broad range of surgical procedures
- Cool LED light eliminates concerns over heat generation
- Vibrant color rendition maximizes visualization of critical anatomy
- 25,000+ hour lamp life limits uncertainty over system brightness and cost of ownership

#### **Greater Freedom**

- Light weight 20' power cord increases freedom from clutter of floor stands and fiber optic cables
- Light weight battery power allows untethered movement from side to side of the operating table and from OR to OR

#### Simple, quiet and comfortable

- Simple 2-point adjustment enables comfortable fit for heads of all sizes
- Unique negative-flow cooling system is whisper-quiet
- Easy to reach 4-position intensity switch is mounted on the headband

# Integra LED Headlight Systems

#### LED Headlight w/Battery & AC/DC Power Supply

Catalog Number	Description
90520US	LED Headlight w/Battery & AC/DC Power Supply - US
90520EU	LED Headlight w/Battery & AC/DC Power Supply - EU
90520UK	LED Headlight w/Battery & AC/DC Power Supply - UK
90520AU	LED Headlight w/Battery & AC/DC Power Supply - AU
90520WUS	LED Headlight w/Battery & AC/DC Power Supply - US, Wide
90520WEU	LED Headlight w/Battery & AC/DC Power Supply - EU, Wide



### LED Headlight w/AC/DC Power Supply Only

Catalog Number	Description
90520DUS	LED Headlight w/AC/DC Power Supply Only - US
90520DEU	LED Headlight w/AC/DC Power Supply Only - EU
90520DUK	LED Headlight w/AC/DC Power Supply Only - UK
90520DAU	LED Headlight w/AC/DC Power Supply Only - AU
90520DWUS	LED Headlight w/AC/DC Power Supply Only - US, Wide
90520DWEU	LED Headlight w/AC/DC Power Supply Only - EU, Wide



# 90520 Integra LED Surgical Headlight System with Batteries and AC/DC Power

Catalog Number	Description	Headlight	Battery Kit	AC/DC Power Kit
90520US	LED Headlight w/Batteries & AC/DC Power Supply - US	90500	90530US	90550US
90520EU	LED Headlight w/Batteries & AC/DC Power Supply - EU	90500	90530EU	90550EU
90520UK	LED Headlight w/Batteries & AC/DC Power Supply - UK	90500	90530UK	90550UK
90520AU	LED Headlight w/Batteries & AC/DC Power Supply - AU	90500	90530AU	90550AU

# 90520D Integra LED Surgical Headlight System with AC/DC Power Only

Catalog Number	Description	Headlight	Battery Kit	AC/DC Power Kit
90520DUS	LED Headlight with AC/DC Power Supply Only - US	90500	N/A	90550US
90520DEU	LED Headlight with AC/DC Power Supply Only - EU	90500	N/A	90550EU
90520DUK	LED Headlight with AC/DC Power Supply Only - UK	90500	N/A	90550UK
90520DAU	LED Headlight with AC/DC Power Supply Only - AU	90500	N/A	90550AU

# **Parts and Accessories**



#### Integra LED Parts and Accessories

Catalog Number	Description
90500	LED Headlight Only
90500W	LED Headband Only, Wide
90500RS	LED Replacement Screw Kit
90542	LED Headband Replacement Pads (10/Box)
90520CS	Nylon LED System Carrying Case
90521	Battery Holster
90522	Lithium Ion Batteries (Qty 2)
90523	Battery Charger, Single Bay
90524	Battery Charger, Dual Bay
90524US	Dual Bay Charger Kit - US
90530US	Battery Kit, US
90530EU	Battery Kit, EU
90530UK	Battery Kit, UK
90530AU	Battery Kit, AU
90550US	AC/DC Power Kit, US
90550EU	AC/DC Power Kit, EU
90550UK	AC/DC Power Kit, UK
90550AU	AC/DC Power Kit, AU
90531	20' DC Power Cord
407421	AC/DC 30W Power Supply
407481	Battery Charger 60W Power Supply (for charger)
90590US	AC Mains Power Cord - US
90590EU	AC Mains Power Cord - EU
90590UK	AC Mains Power Cord - UK
90590AU	AC Mains Power Cord - AU



# MLX 300

### We Are Passionate About Putting Surgery in the Best Light

Sunlight is the ideal light for humans to work. It gives you the color distinction you need in surgery. The MLX Xenon light source gives you true white light, not harsh actinic white light with blue shading.

#### Increased Lamp Warranty - 1000 hours

- True warranty, not pro-rated
- Better cooling leads to longer lamp life

#### **Membrane Touchpad Controls**

- Precise light attenuation in 5% increments
- Lamp and System Age meters

#### **New Floorstand**

- Dedicated to MLX, part number: 001337
- New large casters 75 mm vs. 35 mm
- 10° angle for better viewing of LCD display
- Lower center of gravity

#### **Stand-by Switch**

- System starts up in Stand-by Mode
- Allows electronic self-check
- Put unit into Stand-by to block light when removing the cable
- Starts up at the same light output setting



#### **MLX Light Source And Accessories**

Catalog Number	Description
ooMLX	MLX 300W Xenon light source with 20' power cord
ooMLXEU	MLX 300W Xenon light source with 6.1m EU power cord
ooMLXAU	MLX 300W Xenon light source with 6.1m Australian power cord
ooMLXUK	MLX 300W Xenon light source with 6.1m UK power cord
oo1337	Floorstand for MLX 300W Xenon light source



# UltraLite<sup>®</sup> Pro Headlight

### Maximum Comfort

- Occipital Basket
  - Balances and distributes the weight of the headlight and reduces neck strain
- Cranial Support Stabilizes the headlight
- Adjustable Headband- Provides perfect sizing without temporal compression

### Maximum Headlight Performance:

- Superior optics
- Optimized light output
- Enhanced Design
  - New lighter design
  - 3-Point Support System for comfortable stability
  - Extended linkage for better vertical alignment with surgical loupes

#### Premium Headlight Cable:

- Lightweight and flexible 3mm cable
- Bifurcated design allows even weight distribution and balance
- Translucent outer sheath provides visualization of fiber integrity
- Lifetime Warranty Against Browning

#### **UltraLite Pro Headlight System**

Catalog Number	Description
AX2100BIF	UltraLite Pro headlight with 9ft (275cm) premium bifurcated cable, headlight module and gown clips
001370LX	UltraLite Pro replacement headband
001380LX	Headlight module, black
001388LX9	UltraLite Pro premium bifurcated headlight cable, 9ft (275cm)
001388LX11	UltraLite Pro premium bifurcated headlight cable, 11ft (335cm)
001387LX	UltraLite Pro headlight replacement pad kit



# UltraLite<sup>®</sup> Plus Headlight

### Lightweight ergonomic design

- Our soft, cool UltraLite<sup>®</sup> Plus Sweatband uses a high-tech, lightweight material that's washable and breathable
- Our popular, lightweight UltraLite Plus Headband is ergonomically designed and fully adjustable
- Both designed for maximum comfort during extended wear

#### Superior headlight optics

- Visualize more clearly, see colors more accurately and work more effectively
- Our technology delivers brighter light where you need it: at the center of the spot, not the edges

#### **UltraLite Plus Headlight Systems**

Catalog Number	Description
AX1385BIF	UltraLite Plus Headlight System, sweatband, with extended linkage
AX1375BIF	UltraLite Plus Headlight System, headband, with extended linkage
001290BIF	Replacement Sweatband only
001270BIF	Replacement Headband only
0029009	Replacement Bifurcated Cable only
400282LX	Headband Pad Replacement Kit

#### UltraLite Plus Sweatband



UltraLite Plus Headband



# DLX UltraLite® Pro Camera

#### State of the Art Quality

- The flat beam illumination and custom optics provide superb focus across the entire depth of the visual field
- Super light weight camera adds less than one ounce
- UltraLite Pro design offers 3-point support with occipital basket

#### DLX UltraLite Pro

Catalog Number	Description
003285NTSC	Fixed Focus Headlight Camera System - NTSC (U.S.)
003285PAL	Fixed Focus Headlight Camera System - PAL (outside U.S.)



# **Parts and Accessories**

# Headlight Replacement Components

Catalog Number	Description	Where Used
001270BIF	UltraLite Plus headband only	UltraLite Plus
001290BIF	UltraLite Sweatband only	UltraLite Sweatband
001370LX	UltraLite Pro headband only	UltraLite Pro
001380LX	Mark II Module, black	UltraLite Pro
001328	Mark II Module, gray	UltraLite Pro Module
001190	Mark I Module, gray	UltraLite Plus and Sweatband
001190EL	Mark I Module with extended linkage, gray	UltraLite Plus and Sweatband
005331LX	Light cable clip, teal	UltraLite Plus and Sweatband
001387LX	UltraLite Pro padding	UltraLite Pro
90542	LED padding	LED and Coolvent
004810LX	Gown Clip	All headlights
000286	Joystick	Mark I Module
400282LX	UltraLite Plus padding	UltraLite Plus
400355	Extended linkage kit	All Modules

# Lightsource Replacement Components

Catalog Number	Description	Where Used
001320LX	Complete lamp module	MLX
MLXPSK	Power supply kit	MLX
0025001	Lamp only	9300XSP
000090	Complete lamp module	9300XSP
600972	Power supply kit	9300XSP
001400	Floorstand	9300XSP



Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

• Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.

• Warning: Applicable laws restrict these products to sale by or on the order of a physician.

- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 = 888-980-7742 fax International +1 609-936-5400 = +1 609-750-4259 fax integralife.com/contact

