VENOUS EXTERNAL SUPPORT

← 4.0cm - 23.0cm Ø 3.0 - 8.0mm

Models: 4





Product Description

VEST 2.0 is a cobalt chrome external support intended for the preventive treatment of vein grafts which are being used as conduits in patients who undergo a coronary artery bypass grafting (CABG) procedure.

Product Models

Vein Graft Length	Model by Vein Graft Diameter [PN]	
	ØA	ØB
4.0 cm - 16.0 cm	ØA16 [FG066]	ØB16 [FG069]
16.0 cm - 23.0 cm	ØA23 [FG067]	ØB23 [FG070]

 $3.0 \text{ mm} < \emptyset A \le 5.5 \text{ mm} < \emptyset B \le 8.0 \text{ mm}$

New Features

- VEST 2.0 can be cut to fit the exact length of the vein graft.
- VEST 2.0 is more pliable and more compatible with the vein graft geometry.
- Reduction in number of models. 4 models covering vein grafts lengths of 4.0 - 23.0cm and wide range of diameters between 3.0 - 8.0mm.

Accessories

VEST Selection Tool

A reusable tool for the assessment of vein graft length and diameter, intended for proper selection of the relevant VEST 2.0 model.



Product Code

FG014

VGS Olive Cannula

A reusable cannula for easier implantation of the VEST 2.0 device over the vein graft.



Product Code

FG013



