


EXTERNAL SUPPORT FOR PERIPHERAL VASCULAR RECONSTRUCTION

→ 30cm - 90cm

Ø 3.5 - 8.0mm

Models: 12





Product Description

FRAME is a kink resistant, cobalt-chrom, external support for vein grafts in peripheral bypass and reconstruction procedures, which mitigates the pathological remodeling of vein grafts.

Product Models

Vein Graft External Diameter		30cm	60cm	90cm
ØA	3.5 - 4.5mm	FG031	FG034	FG037
ØB	4.6 - 5.5mm	FG040	FG043	FG046
ØC	5.6 - 6.5mm	FG049	FG052	FG055
ØD	6.6 - 8.0mm	FG058	FG061	FG064



— Can be cut to perfectly fit the graft length

Accessories

FRAME Selection Tool

A reusable tool for the assessment of vein graft diameter, intended for proper selection of the relevant FRAME model.



Product Number

FG073

VGS Olive Cannula

A reusable cannula for easier implantation of the FRAME device over the vein graft.



Product Number

FG013



