# FRAME FR<sup>TM</sup>

EXTERNAL SUPPORT FOR HIGH-FLOW & ANERUSYMAL AV FISTULA RECONSTRUCTION

Up to 20cm

Ø 4.6 - 8.0mm

Models: 6



## **Product Description**

FRAME FR introduces a new approach for repair and stabilization of high flow and aneurysmal AVF. Deployment of FRAME FR over a reconstructed venous segment reduces and maintains AVF flow within desired range, and prevents re-dilatation.

## FRAME FR Models

FRAME FR Model	Vein External Diameter	10cm	20cm
ØB	4.6 - 5.5mm	FG078	FG081
ØC	5.6 - 6.5mm	FG079	FG082
ØD	6.6 - 8.0mm	FG080	FG083



Can be cut to perfectly fit the graft length

### Accessories

#### FRAME Selection Tool

A reusable tool for the assessment of vein graft diameter, intended for proper selection of the relevant FRAME FR model.



**Product Number** 

#### VGS Olive Cannula

A reusable cannula for easier implantation of the FRAME FR device over the vein graft.



**Product Number** 

FG013

#### FRAME FR Mandrel Set

A set of three reusable mandrels for the reconstruction of the vein segment while maintaining effective lumen. The mandrels are compatible with the FRAME FR models.





MedNet EC-REP GmbH Borkstrasse 10, 48163 Münster, Germany



