

Swartz™ Braided Transseptal Guiding Introducers

SL Series

63 cm Length

8 F–8.5 F


Product Highlights

- Braided sheath facilitates torqueability, pushability and catheter stability
- Ultra-soft tip for atraumatic performance
- Multiple curve options available
- Smaller offset sideholes reduce likelihood of guidewire exiting sidehole
- Ultimium™ valve for effective hemostasis and minimal risk of air inspiration
- Accommodates a 71 cm BRK™ Transseptal Needle (Reorder Number 407200 or 407201)

Ordering Information

Contents: Swartz™ Braided SL Transseptal Guiding Introducer Sheath with Hemostasis Valve and Sideport, Dilator, and 180 cm Super Stiff Guidewire with Finger Straightenable 3 mm “J” (1 unit per box)

Reorder Number	Sheath French Size	Dilator French Size	Maximum Guidewire Diameter (in)	Curve Type	Sheath Usable Length (cm)	Dilator Usable Length (cm)
407449	8	8	.032	SL0™	63	67
407439	8	8	.032	SL1™	63	67
407441	8	8	.032	SL2™	63	67
407443	8	8	.032	SL3™	63	67
407446	8	8	.032	SL4™	63	67
407451	8.5	8.5	.032	SL0™	63	67
407453	8.5	8.5	.032	SL1™	63	67
407455	8.5	8.5	.032	SL2™	63 <td 67	
407457	8.5	8.5	.032	SL3™	63	67
407459	8.5	8.5	.032	SL4™	63	67



Swartz™ Braided Transseptal Guiding Introducers

SL Series

81 cm Length

8 F–8.5 F


Product Highlights

- Braided sheath facilitates torqueability, pushability and catheter stability
- Ultra-soft tip for atraumatic performance
- Multiple curve options available
- Smaller offset sideholes reduce likelihood of guidewire exiting sidehole
- Ultimum™ valve for effective hemostasis and minimal risk of air inspiration
- Accommodates an 89 cm BRK™ Transseptal Needle (Reorder Number 407205)

Ordering Information

Contents: Swartz™ Braided SL Transseptal Guiding Introducer Sheath with Hemostasis Valve and Sideport, Dilator, and 180 cm Super Stiff Guidewire with Finger Straightenable 3 mm “J” (1 unit per box)

Reorder Number	Sheath French Size	Dilator French Size	Maximum Guidewire Diameter (in)	Curve Type	Sheath Usable Length (cm)	Dilator Usable Length (cm)
407450	8	8	.032	SL0™	81	85
407440	8	8	.032	SL1™	81	85
407442	8	8	.032	SL2™	81	85
407452	8.5	8.5	.032	SL0™	81	85
407454	8.5	8.5	.032	SL1™	81	85
407456	8.5	8.5	.032	SL2™	81	85



Information contained herein for distribution outside of the U.S. only.

Electrophysiology International Catalogue

©2013 St. Jude Medical. All Rights Reserved.

Swartz™ Braided Transseptal Guiding Introducers LAMP™ Series 63 cm Length 8 F–8.5 F


Product Highlights

- Braided sheath facilitates torquability, pushability and catheter stability
- Ultra-soft tip for atraumatic performance
- Three curve options available (LAMP™ 45, LAMP™ 90, LAMP™ 135)
(LAMP = Left Atrial Multi Purpose)
- Smaller offset sideholes reduce likelihood of guidewire exiting sidehole
- Ultimum™ valve for effective hemostasis and minimal risk of air inspiration
- Accommodates a 71 cm BRK™ Transseptal Needle
(Reorder Number 407205)

Ordering Information

Contents: Swartz™ Braided LAMP™ Transseptal Guiding Introducer Sheath with Hemostasis valve and Sideport, Dilator, and 180 cm Super Stiff Guidewire 3 mm “J” (1 unit per box)

Reorder Number	Sheath French Size	Dilator French Size	Maximum Guidewire Diameter (in)	Curve Type	Sheath Usable Length (cm)	Dilator Usable Length (cm)
407360	8	8	0.032	LAMP™ 45	63	67
407362	8.5	8.5	0.032	LAMP™ 45	63	67
407356	8	8	0.032	LAMP™ 90	63	67
407358	8.5	8.5	0.032	LAMP™ 90	63	67
407364	8	8	0.032	LAMP™ 135	63	67
407366	8.5	8.5	0.032	LAMP™ 135	63	67



Swartz™ Braided Transseptal Guiding Introducers

LAMP™ Series

81 cm Length

8 F–8.5 F


Product Highlights

- Braided sheath facilitates torquability, pushability and catheter stability
- Ultra-soft tip for atraumatic performance
- Three curve options available (LAMP™ 45, LAMP™ 90, LAMP™ 135)
(LAMP = Left Atrial Multi Purpose)
- Smaller offset sideholes reduce likelihood of guidewire exiting sidehole
- Ultim™ valve for effective hemostasis and minimal risk of air inspiration
- Accommodates an 89 cm BRK™ Transseptal Needle
(Reorder Number 407205)

Ordering Information

Contents: Swartz™ Braided LAMP™ Transseptal Guiding Introducer Sheath with Hemostasis valve and Sideport, Dilator, and 180 cm Super Stiff Guidewire 3 mm “J” (1 unit per box)

Reorder Number	Sheath French Size	Dilator French Size	Maximum Guidewire Diameter (in)	Curve Type	Sheath Usable Length (cm)	Dilator Usable Length (cm)
407363	8.5	8.5	0.032	LAMP™ 45	81	85
407357	8	8	0.032	LAMP™ 90	81	85
407359	8.5	8.5	0.032	LAMP™ 90	81	85
407365	8	8	0.032	LAMP™ 135	81	85
407367	8.5	8.5	0.032	LAMP™ 135	81	85



Information contained herein for distribution outside of the U.S. only.

Electrophysiology International Catalogue

©2013 St. Jude Medical. All Rights Reserved.