

Pacel™ Bipolar Pacing Catheter,  
Ventricular Pacing 60 Degree Curve  
Torque Control

Product Highlights

- Designed to provide reliable pacing and sensing in the ventricle from the superior vena cava approach
- Precise control without excessive stiffness
- Optimal torque control for challenging anatomy
- Solid core enables user to customize the curve
- 10 cm depth markers to enhance placement
- Smooth, hemocompatible polyurethane surface
- Stainless steel center conductor and electrodes
- Clearly marked terminals
- Shrouded 2 mm pin connector leads

Ordering Information

Shrouded Pin Contents: Pacing Catheter and Pin Adapters (5 units per box)

Exposed Pin Contents: Pacing Catheter (5 units per box)

Kit Contents: Pacing Catheter and Needle Peel-Away Introducer (5 units per box)

Shrouded Pin Reorder Number	Exposed Pin Reorder Number**	French Size	Electrode Spacing (cm)	Usable Length (cm)	Recommended Introducer Size (F)
401765	401664	4	1	110	4
401775	401686	5	2.5	110	5
401766	401665	5	1	110	5
401776	401688	6	2.5	110	6
401767	401666	6	1	110	6
401768	401672	7	1	110	7
N/A	401641 (Kit)	5	1	110	5
N/A	401645 (Kit)	5	2.5	110	5
N/A	401642 (Kit)	6	1	110	6
N/A	401646 (Kit)	6	2.5	110	6



\*Allow up to 3 weeks for delivery.

\*\*Not available in countries requiring CE Mark.

Product referenced is approved for CE Mark.

# Pacel™ Bipolar Pacing Catheter, Ventricular Pacing with Right Heart Curve Torque Control

## Product Highlights

- Designed to provide reliable pacing and sensing in the ventricle from the inferior vena cava approach
- Precise control without excessive stiffness
- Optimal torque control for challenging anatomy
- Solid core enables user to customize the curve
- 10 cm depth markers to enhance placement
- Smooth, hemocompatible polyurethane surface
- Stainless steel center conductor and electrodes
- Clearly marked terminals
- Shrouded 2 mm pin connector leads

## Ordering Information

Shrouded Pin Contents: Pacing Catheter and Pin Adapters (5 units per box)

Exposed Pin Contents: Pacing Catheter (5 units per box)

Kit Contents: Pacing Catheter and Needle Peel-Away Introducer (5 units per box)

Shrouded Pin Reorder Number	Exposed Pin Reorder Number**	French Size	Electrode Spacing (cm)	Usable Length (cm)	Recommended Introducer Size (F)
401769	401673	4	1	110	4
401770	401674	5	1	110	5
401771	401675	6	1	110	6
401772	401676	7	1	110	7
N/A	401643 (Kit)	5	1	110	5
N/A	401644 (Kit)	6	1	110	6



\*Allow up to 3 weeks for delivery.

\*\*Not available in countries requiring CE Mark.

Product referenced is approved for CE Mark.

## Pacel™ Bipolar Pacing Catheter, Atrial “J” Pacing Torque Control

### Product Highlights

- Designed to provide reliable pacing and sensing in the right atrium from the superior vena cava approach
- Precise control without excessive stiffness
- Optimal torque control for challenging anatomy
- Solid core enables user to customize the curve
- 10 cm depth markers to enhance placement
- Smooth, hemocompatible polyurethane surface
- Stainless steel center conductor and electrodes
- Clearly marked terminals
- Shrouded 2 mm pin connector leads

### Ordering Information

Shrouded Pin Contents: Pacing Catheter and Pin Adapters (5 units per box)

Exposed Pin Contents: Pacing Catheter (5 units per box)

Shrouded Pin Reorder Number	Exposed Pin Reorder Number**	French Size	Electrode Spacing (cm)	Usable Length (cm)	Recommended Introducer Size (F)
401773	401680	5	1	110	5
401774	401682	6	1	110	6



\*Allow up to 3 weeks for delivery.

\*\*Not available in countries requiring CE Mark.

Product referenced is approved for CE Mark.

## Pacel™ Bipolar Pacing Catheter, Straight Flow Directed

### Product Highlights

- Designed to provide quick, reliable temporary pacing in emergency situation with or without the use of fluoroscopy
- Flexible design for delicate placement of catheter
- Rapid inflation and deflation of balloon
- 10 cm depth markers to enhance placement
- Smooth, hemocompatible polyurethane surface
- Clearly marked terminals
- Shrouded 2 mm pin connector leads

### Ordering Information

401761 Contents: Pacing Catheter, Syringe and Pin Adapters (1 unit per box)

401615 Contents: Pacing Catheter and Syringe (1 unit per box)

401762 Kit Contents: Pacing Catheter, Syringe, Pin Adapters and Needle Peel-Away Introducer (1 unit per box)

401616 Kit Contents: Pacing Catheter, Syringe, Needle Peel-Away Introducer and ECG Adapter (1 unit per box)

Shrouded Pin Reorder Number	Exposed Pin Reorder Number*	French Size	Electrode Spacing (cm)	Usable Length (cm)	Recommended Introducer Size (F)
401761	401615*	5	1	110	6
401762 (Kit)	401616 (Kit)	5	1	110	6



\*Allow up to 3 weeks for delivery.

\*Not available in countries requiring CE Mark.

Product referenced is approved for CE Mark.

# Pacel™ Bipolar Pacing Catheter, Right Heart Curve

### Flow Directed

## Product Highlights

- Designed to provide quick, reliable temporary pacing in emergency situation with or without the use of fluoroscopy
- Flexible design for delicate placement of catheter
- Rapid inflation and deflation of balloon
- 10 cm depth markers to enhance placement
- Smooth, hemocompatible polyurethane surface
- Clearly marked terminals
- Shrouded 2 mm pin connector leads

## Ordering Information

401763 Contents: Pacing Catheter, Syringe and Pin Adapters (1 unit per box)

401623 Contents: Pacing Catheter and Syringe (1 unit per box)

401764 Kit Contents: Pacing Catheter, Syringe, Pin Adapters and Needle Peel-Away Introducer (1 unit per box)

401624 Kit Contents: Pacing Catheter, Syringe, Needle Peel-Away Introducer and ECG Adapter (1 unit per box)

Shrouded Pin Reorder Number	Exposed Pin Reorder Number*	French Size	Electrode Spacing (cm)	Usable Length (cm)	Recommended Introducer Size (F)
401763	401623*	5	1	110	6
401764 (Kit)	401624 (Kit)	5	1	110	6



\*Allow up to 3 weeks for delivery.

\*Not available in countries requiring CE Mark.

Product referenced is approved for CE Mark.