

Merlin™

Patient Care System (PCS)

Product Highlights

- 15" touch screen clearly displays programming and diagnostic screens
- New user interface allows for faster patient management
- Continuous, simultaneous display of surface ECGs, intracardiac electrograms (EGMs) and annotated event markers allow quick interpretation
- Built-in top-load/top-exit printer quickly and quietly prints full-page (8-inch) reports for patient charts
- Integrated cable storage speeds setup and saves space with always-connected cables and ample storage space




Ordering Information

Model Number	Part Number	Description
3650		Merlin Patient Care System (PCS)
Merlin PCS accessories		
Model Number	Part Number	Description
3001		3-Lead ECG Patient Cable
3626		5-Lead ECG Patient Cable
3134	60000909-001	VGA Cable and Adapter (female to male) (25' length)
3615	60004294-001	Adapter for 3150 PSA Wand (required for use of PSA Wand Model 3150 with the Merlin Patient Care System)
3616	60005260-001	Wand Extension Cable (4' length)
3617	60005251-001	External ECG Input Cable (25' length)
3620	60005254-001	External Floppy Disk Drive
EX3621-2GB	100006806	Flash Drive (2 GB)
3622	60005256-001	Shoulder Strap
3623	60005257-001	USB to RS232 Serial Adapter (for direct Paceart™ connection)
3630	60002876-001	Merlin Patient Care System Telemetry Wand (with Magnet)
3630M	60002876-097	Magnet
3638	50019403-001	Antenna (required to enable radio frequency [RF] communication between Merlin PCS and St. Jude Medical™ implantable devices with RF communication capability)
3643	60003605-001	Thermal Paper

Customer Support: 46-8-474-4756

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use. Devices depicted may not be available in all countries. Check with your St. Jude Medical representative for product availability in your country. Unless otherwise noted, ™ indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2011 St. Jude Medical, Inc. All rights reserved.



Customer Support: 46-8-474-4756

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use. Devices depicted may not be available in all countries. Check with your St. Jude Medical representative for product availability in your country. Unless otherwise noted, TM indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2011 St. Jude Medical, Inc. All rights reserved.

Item GMC827EN



ST. JUDE MEDICAL
MORE CONTROL. LESS RISK.