Attune™ Flexible Adjustable Annuloplasty Ring

Product Highlights


- · Fully flexible design allows dynamic movement of the annulus
- Symmetrical and asymmetrical adjustability is designed to achieve ideal leaflet coaptation and match each patient's annulus
- Adjustability after placement is designed to help eliminate residual mitral regurgitation
- Independent adjustability allows the annuloplasty to be localized to one side
- Suitable for implant using open sternotomy, minimally invasive and robotic approaches

Ordering Information


Contents: Flexible Adjustable Annuloplasty Ring (1 unit per box)

Reorder Number	Ring Size	Fixed Dimension (mm)	Ring Inner Dimension (mm)	Ring Outer Dimension (mm)
AFR-25	25	25	28.4	33.4
AFR-43	43	43	51.6	56.6
AFR-27	27	27	31.2	36.2
AFR-29	29	29	34.0	39.0
AFR-31	31	31	36.6	41.6
AFR-33	33	33	39.4	44.4
AFR-35	35	35	40.9	45.9
AFR-37	37	37	43.7	48.7
AFR-39	39	39	46.2	51.2
AFR-41	41	41	49.0	54.0





A: Fixed Dimension (mm) B: Ring Inner Dimension (mm) C: Ring Outer Dimension (mm)



Rx Only

Please review the Instructions for Use prior to using these devices for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

Unless otherwise noted, [®] or [™] indicates a registered or unregistered trademark or service mark owned by, or licensed to, St. Jude Medical, Inc. or one of its subsidiaries. ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies.

©2013 St. Jude Medical, Inc. All rights reserved.

100033292 - May 2012

Customer Support: 855-478-5833

Last Updated: April 19, 2013