Pericardial Patch WITH ENCAPTM AC TECHNOLOGY

ENHANCED BIOCOMPATIBILITY FOR LASTING PERFORMANCE

The Pericardial Patch with EnCap AC Technology* combines strong, durable bovine pericardium with a proprietary anti-calcification treatment, making it suitable for a variety of cardiac repairs while offering improved handling and enhanced biocompatibility.



IMPROVED HANDLING AND SUTURABILITY SUPPORT CARDIOVASCULAR REPAIR

- Ready-to-use, rinseless preparation saves time during procedures.
- Bovine pericardium provides improved handling and suturability compared with synthetic patches.¹
- The strength of glutaraldehyde-fixed tissue enhances durability and helps resist undesirable changes such as patch shrinkage and aneurysm formation, even in high-stress repairs.²⁻⁶
- Soft, pliable tissue conforms to anatomy and sutures into place with minimal leaking along suture line.

ANTI-CALCIFICATION TREATMENT ENHANCES BIOCOMPATIBILITY AND DURABILITY

- Proprietary EnCap AC Technology caps residual aldehydes to reduce antigenicity and cytotoxicity.^{57,8}
- Resists calcification and promotes rapid binds, thorough healing with endothelial cell covering.⁷⁻¹⁰
- Improved endothelialization strengthens the reconstruction or repair, helping reduce calcification and other degeneration.⁷⁻⁹

APPROPRIATE FOR A WIDE RANGE OF CARDIAC AND VASCULAR REPAIRS⁸

- Annular reconstruction³
- Endocarditis leaflet repairs
- Septal defect repairs
- Aortic root enlargement
- Other vascular repairs.



ORDERING INFORMATION

Pericardial Patch

Model Number	Patch Size (cm)	Nominal Thickness (mm)
C0205	2 x 5	0.20 - 0.40
C0405	4 x 5	0.15 - 0.25
C0510	5 x 10	0.20 - 0.40
C0914	9 x 14	0.20 - 0.40

All sizes not currently available in all markets.

References:

- Crawford FA Jr, Sade RM, Spinale F. Bovine pericardium for correction of congenital heart defects. Ann Thorac Surg. 1986;41(6):602-5. 1.
- 2.
- 3. 4
- Frater RWM, Vetter HO, Zussa C, et al. Chordal replacement in mitral valve repair. Circulation. 1990;82[suppl IV]:IV-125-IV-130. David TE, Feindel CM, Armstrong, S, et al. Reconstruction of the mitral annulus: a ten-year experience. J Thorac Cardiovasc Surg. 1995;110(5):1323-32. Bjornstad K, Duran RM, Nassau KG, et al. Clinical and echocardiographic follow-up after aortic valve reconstruction with bovine or autologous pericardium. Am Heart J. 1996;132(6):1173-8. Gong G, Seiffer E, WD Lyman, et al. Bioprosthetic cardiac valve degeneration: role of inflammatory and immune reactions. J Heart Valve Dis. 1993;2(6):684-93. 5.
- Gong G, Ling Z, Seifter E, et al. Aldehyde tanning: the villain in bioprosthetic calcification. Eur J Cardiothorac Surg. 1991;5:288-99. Frater RWM, Seifter E, Liao K, et al. Anticalcification, proendothelial, and anti-inflammatory effect of post-aldehyde polyol treatment of bioprosthetic material. In: Gabbay S, Wheatley DJ (eds.). 6. 7.
- Advances in Anticalcific and Antidegenerative Treatment of Heart Valve Bioprostheses. Austin, TX: Silent Partners Inc; 1997:105-14. Frater RWM, Liao K, Seifter E. Stentless chordally supported mitral bioprosthetic valve. In: Gabbay S, Frater RWM (eds.) New Horizons and the Future of Heart Valve Bioprostheses. Austin, TX: 8. Silent Partners Inc; 1994:103-19.
- Hoffman D, Gong G, Liao K, et al. Spontaneous host endothelial growth on bioprostheses. Circulation. 1992;86[suppl II]:II-75-II-79.
- 10 Moritz A, Grimm M, Eybl E, et al. Improved spontaneous endothelialization by postfixation treatment of bovine pericardium. Eur J Cardiothorac Surg. 1991;5:155-9.

Abbott Vascular International BVBA



Park Lane, Culliganlaan 2b, 1831 Diegem, Belgium Products intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use provided inside the product carton (when available), at eifu abbottvascular.com or at manuals.sjm.com for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. Photo(s) on file at Abbott. Information contained herein is for distribution for Europe, Middle East and Africa ONLY. Please check the regulatory status of the device before distribution in areas where CE marking is not the regulation in force. For more information, visit our website at www.abbott.com © 2018 Abbott. All Rights Reserved. 9-EH-I-8616-01 10-2018

