

		Off Time (minutes)								
		0.2	0.3	0.5	0.8	1.1	1.8	3	5	10
On Time (seconds)	7	58%	44%	30%	20%	15%	10%	6%	4%	2%
	14	69	56	41	29	23	15	9	6	3
	21	76	64	49	36	29	19	12	8	4
	30	81	71	57	44	35	25	16	10	5
	60	89	82	71	59	51	38	27	18	10

- Recommended
- Not recommended*

- Recommended progression for duty cycle
- Not available if Cardiac Based Seizure Detection (CBSD) feature is enabled