## PressureWire<sup>™</sup> Aeris<sup>™</sup>

## FFR Measurement System

## **Product Highlights**

- Wirelessly measures FFR and is designed to offer exceptional control and maneuverability even in tortuous vessels
- Used with PressureWire Receiver (sold seperately)
- Together with a software upgrade to the cath lab hemodynamic recording system it forms a fully integrated FFR measurement system
- Patented sensor technology enables accurate FFR measurement
- A hydrophilic polymer sleeve isolates the wires cabling and electronics to increase signal reliability

## Ordering Information

Contents: PressureWire Aeris (5 units per box)

Reorder		
Number	Length (cm)	
12056	PressureWire Aeris 175	



Product referenced is approved for CE Mark.

