

# PressureWire™ Certus™ FFR Measurement System

## Product Highlights

- A continuous taper core wire provides excellent torque response for excellent maneuverability in tortuous anatomies
- Used together with RadiAnalyzer™ Xpress (sold separately)
- Optimized sensor-to-wire transition enhances ability to navigate sharp-angled side branches and makes less sensitive to prolapse
- Patented high fidelity sensor technology for simultaneous measurement of pressure (FFR), flow (CFR), microcirculatory resistance (IMR\*) and intravascular temperature
- Together with RadiAnalyzer™ Xpress forms an advanced measurement system
- Polymer sleeve isolates cabling to protect electronics increasing signal reliability

## Ordering Information

Contents: PressureWire Certus (5 units per box)

Reorder Number	Length (cm)
12006	PressureWire Certus 175



Product referenced is approved for CE Mark.