

EnligHTN™

Renal Artery Ablation Catheter



THE POWER OF A PREDICTABLE PATTERN

The EnligHTN renal denervation system is the first, multi-electrode ablation technology to deliver a pre-determined ablation set in a single positioning. Four electrodes on the EnligHTN ablation catheter administer a repeatable pattern of evenly spaced ablations with each placement. By reducing the number of catheter manipulations needed to disrupt the sympathetic nerve network, this new technology offers cutting-edge efficiencies that enable unparalleled speed and placement accuracy in renal denervation.¹⁻³

THE BENEFITS OF MULTI-ELECTRODE ABLATION

FAST

- The EnligHTN ablation catheter administers a consistent pattern of four ablations with each placement — in 90-second intervals — potentially reducing renal denervation procedure time.^{2,3}

EFFICIENT

- Unique multi-ablation basket design provides four ablations without readjustment.

PRECISE

- Consistent, four-point contact with arterial wall designed to increase procedural accuracy.

MINIMIZED EXPOSURE

- Minimal catheter repositioning may result in shorter procedure time, reduced contrast and fluoroscopic exposure.

ENLIGHTN MULTI-ELECTRODE ABLATION¹



Precision Pattern

SINGLE-ELECTRODE ABLATION¹



Random Pattern



ST. JUDE MEDICAL™

MORE CONTROL. LESS RISK.

YOUR PARTNER IN RENAL DENERVATION TREATMENT

The EnligHTN renal denervation system offers advanced renal artery ablation treatment and is designed to deliver radiofrequency (RF) energy to the renal nerves to achieve targeted denervation. The core system components are the EnligHTN RF ablation generator, the EnligHTN renal artery ablation catheter and the EnligHTN guiding catheter.

Ordering Information

EnligHTN Renal Denervation System	
Reorder Number	Product Name
ENL-LG-01	EnligHTN Renal Artery Ablation Catheter - large
ENL-SM-01	EnligHTN Renal Artery Ablation Catheter - small
ENL-8F55	EnligHTN Guiding Catheter
ENL-T115	EnligHTN RF Ablation Generator
ENL-CC-01	EnligHTN 1641 Connector Cable
PREP-01	Patient Return Electrode
ENL-PC-01	EnligHTN Power Cable - replacement



EnligHTN Renal Artery Ablation Catheter



EnligHTN Guiding Catheter



EnligHTN RF Ablation Generator

1. St. Jude Medical. Data on file.
2. St. Jude Medical Instruction for Use, Ablation catheter #10006654702.
3. Medtronic Symplicity Instruction for Use.

ATRIAL FIBRILLATION

CARDIAC RHYTHM MANAGEMENT

CARDIOVASCULAR

NEUROMODULATION

Global Headquarters

One St. Jude Medical Drive
St. Paul, Minnesota 55117
USA
+1 651 756 2000
+1 651 756 3301 Fax

Cardiovascular Division

177 East County Road B
St. Paul, Minnesota 55117
USA
+1 651 756 4470
+1 651 756 4466 Fax

SJM Coordination Center BVBA

The Corporate Village
Da Vincilaan 11 Box F1
1935 Zaventem, Belgium
+32 2 774 68 11
+32 2 772 83 84 Fax

St. Jude Medical (Hong Kong) Ltd.

Suite 1608, 16/F Exchange Tower
33 Wang Chiu Road
Kowloon Bay, Kowloon
Hong Kong SAR
+852 2996 7688
+852 2956 0622 Fax

St. Jude Medical Brasil Ltda.

Rua Itapeva, 538
5º ao 8º andar
01332-000 – São Paulo – SP
Brazil
+55 11 5080 5400
+55 11 5080 5423 Fax

St. Jude Medical Japan Co., Ltd.

Shiodome City Center 15F
1-5-2, Higashi-Shinbashi
Minato-ku
Tokyo 105-7115
Japan
+81 3 6255 6370
+81 3 6255 6371 Fax

SJMenligHTN.com
SJMprofessional.com



ST. JUDE MEDICAL™
MORE CONTROL. LESS RISK.

Product referenced is approved for CE Mark.

EnligHTN, ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2012 St. Jude Medical, Inc. All rights reserved.

Item GMVASC400EN IPN 2155-12

