

# IBI-1500T11 Cardiac Ablation Generator


## Product Highlights

- High-power output up to 150 watts
- Temperature and power controlled generator
- Two way communication with the Cool Point™ Irrigation Pump

## Ordering Information

Contents: 150-watt RF generator (AC). Includes: RF generator, generator power cord, 1641 connecting cable, 1804-S connecting cable, 1452 footswitch, 1710 grounding cable, electrosurgical patient plates (2 each), Data Logger software and Operator's Manual

Reorder Number	Model Number	Description
85523	IBI-1500T11	150-watt RF generator



# Cool Point™ Irrigation Pump

## Product Highlights

- Peristaltic pump
- Two-way communication with the IBI-1500T11 Cardiac Ablation Generator
- 2 microliter bubble detectors
- In-line pressure sensor to detect downstream occlusion
- Compact, lightweight design
- Ergonomic carrying handle
- Includes: Irrigation Pump, Power Cord, Tubing Set, Pole Clamp, 1779 Connecting Cable, Operator’s Manual

## Ordering Information


Reorder Number	Description
85784	Cool Point Irrigation Pump

### Required Cool Point Tubing Set

Reorder Number	Description	Length (cm)
85785	Cool Point Tubing Set (sold individually)	260

### Required Cool Point Accessories

Reorder Number	Description	Length (cm)
85786	Communications cable, 20-feet, between Cool Point Irrigation Pump and IBI-1500T11 Cardiac Ablation Generator	260



## IBI-1500T Series Cardiac Ablation Generator System Accessories

### Product Highlights

- Compatible cables and accessories for IBI-1500T Series Cardiac Ablation Generator

### Ordering Information

Reorder Number	Model Number	Description
85452	1452	Foot switch with cover and cable, (2.4m)
85453	1453	Foot switch with cover and cable (5m)
85454	1457	Foot switch with cover and cable (10m)
85455	1455	Foot switch with cover and cable (15m)
85524	N/A	IBI-1500T11 Remote Control
85526	N/A	IBI-1500T11 Remote Control Cable (6m)
85527	N/A	IBI-1500T11 Remote Control Cable (10m)
85528	N/A	IBI-1500T11 Remote Control Cable (15m)
85641	1641	Connecting cable from IBI-1500T Series Cardiac Ablation Generator to ablation catheter
85710	1710	250 cm grounding cable
85731	1731	Electrosurgical patient plate
85809	1804-S	Connecting cable

## Irrigated RF Therapy - Accessories

### Required Cool Point Tubing Set

Reorder Number	Description	Length (cm)
85785	Cool Point™ Tubing Set (sold individually)	260

### For Use With Therapy Cool Path Ablation Catheters and Stockert Equipment

Reorder Number	Description	Length (cm)
85787	Irrigation Tubing Set	319