

Therapy™ Cool Flex™ Ablation Catheter

Irrigated 4 mm Tip Thermocouple

Quadripolar

7 F

Product Highlights

- Flexible tip electrode can conform to tissue
- Laser cut tip designed for optimal irrigation flow profile
- 0.5 mm spacing between tip and first band electrode to enhance signal quality
- Uni-directional steering
- Push/pull handle

Ordering Information

Contents: 7 F irrigated ablation catheter (1 unit per box)

| Reorder Number | Description | Electrode Spacing (mm) | Tip Electrode (mm) | Band Electrode (mm) | Curve Type | Usable Length (cm) |
|----------------|-------------------------------|------------------------|--------------------|---------------------|------------|--------------------|
| A088015 | 1304-CF-7-0.5(5)2-M-TE4BE1EB | 0.5-5-2 | 4 | 1 | Medium | 110 |
| A088016 | 1304-CF-7-0.5(5)2-L-TE4BE1EB | 0.5-5-2 | 4 | 1 | Large | 110 |
| A088017 | 1304-CF-7-0.5(5)2-XL-TE4BE1EB | 0.5-5-2 | 4 | 1 | X-Large | 110 |
| A088018 | 1304-CF-7-0.5(5)2-L1-TE4BE1EB | 0.5-5-2 | 4 | 1 | L1 | 110 |


Required Catheter Connecting Cable – Page 165

Reorder Number

SJM: 85641

Stockert: 85713, 85709

Medtronic: 85711, 85708



Information contained herein for distribution outside of the U.S. only.

Electrophysiology International Catalogue

©2013 St. Jude Medical. All Rights Reserved.

