

# Inquiry™ Cardioversion Electrophysiology Catheter

## Single Catheter System

### 7F Steerable Curve


### Ordering Information

Contents: 7F 14-Pole cardioversion catheter (1 unit per box)

Item Number	Description	Tip (mm)	Band (mm)	Spacing (mm)	Curve Type
81998	1114-CV-7-3(80)3-L-TE4BE4	4	4	3-80-3	Large

### Required Catheter Connecting Cables

Item Number	Model Number	Description	Length (m)
85201	2101-C	Connecting cable for cardioversion catheter, single catheter system	2.4



# Inquiry™ Cardioversion Electrophysiology Catheter Dual Catheter System 6F Fixed Curve

## Ordering Information

Contents: **(1 unit per box)**

Item Number	Description	Tip (mm)	Band (mm)	Spacing (mm)	Curve Type
80993	1010-CV-6-5-J-TE4BE4	4	4	5	J
80994	1014-CV-6-2-J-TE4BE4	4	4	2	J

## Required Catheter Connecting Cables

Item Number	Model Number	Description	Length (m)
85205	2102-C	Connecting cable for cardioversion catheter, dual catheter system	2.4



# Inquiry™ Cardioversion Electrophysiology Catheter Dual Catheter System 6F Steerable Curve

## Ordering Information

Contents: 6F 10-Pole cardioversion catheter (1 unit per box)

Item Number	Description	Tip (mm)	Band (mm)	Spacing (mm)	Curve Type
81996	1110-CV-6-2-XL-TE4BE4	4	4	2	X-Large
81997	1110-CV-6-2-L-TE4BE4	4	4	2	Large

## Required Catheter Connecting Cables

Item Number	Model Number	Description	Length (m)
85205	2102-C	Connecting cable for cardioversion catheter, dual catheter system	2.4

