

# Intracardiac Imaging ViewFlex™ PLUS ICE Catheter

## Description

The ViewFlex PLUS ICE Catheter is a temporary intracardiac ultrasound catheter intended for use in patients to visualize cardiac structures and blood flow within the heart when connected to a ViewMate™ ultrasound console via the ViewFlex Catheter Interface Module.

## Product Highlights

- 64-element phased array transducer
- Single-handed deflection mechanism
- Large curvature radius
- Steering angles to 120 degrees
- Reliable tip stability

## Ordering Information

Reorder Number	Description
100026714	ViewFlex PLUS ICE Catheter for ViewMate™ Z Intracardiac Ultrasound Console
100042824	ViewFlex PLUS ICE Catheter for ViewMate™ II System
100043720	ViewFlex PLUS ICE Catheter Interface Module for ViewMate Z Intracardiac Ultrasound Console
100038191	ViewFlex PLUS ICE Catheter Interface Module for ViewMate II System

