

ViewFlex™ Xtra Intracardiac Echocardiography Catheter



STREAMLINE PROCEDURES SINGLE-HANDEDLY

The ViewFlex Xtra ICE Catheter provides seamless maneuverability, allowing physicians to focus more on the procedure and less on catheter manipulation to enhance procedural workflow.

Seamless Control and Maneuverability


- Ergonomic handle with auto-lock steering allows for single-handed operation, which can minimize catheter repositioning
- Four-way tip deflection enables optimal viewing planes
- Extended catheter reach on the distal curve assists in right ventricular access

Exceptional Intracardiac Imaging

- The combination of the ViewFlex Xtra catheter and ViewMate™ Z Ultrasound Console provides crisp, clear images
- Vivid color Doppler display provides accurate interpretation of blood flow direction and velocity

Enhanced Lab Workflow

- An integrated catheter connecting cable eliminates the inconvenient step of placing a sterile sleeve during procedure setup
- The lightweight and compact footprint of the ViewMate Z Console requires less space for improved workflow in the lab
- The catheter interface module optimizes placement of the ultrasound system in the lab, reducing obstructions near the patient table



ViewFlex Xtra ICE Catheter Specifications

Feature	Specification
Catheter size	9 F
Introducer size	10 F
Insertable length	90 cm
Tip Deflection	Anterior/Posterior, Left/Right +/- 120 degree flexion
Transducer type	64 element phased array
Viewing sector	80 degrees

Ordering Information

Reorder Number	Product Description
100046962	ViewFlex Xtra ICE Catheter - Structural Heart
100046963	ViewFlex Xtra ICE Catheter - Electrophysiology
100043720	ViewFlex Catheter Interface Module for ViewMate Z Console
100038191	ViewFlex Catheter Interface Module for ViewMate II System

ATRIAL FIBRILLATION CARDIAC RHYTHM MANAGEMENT CARDIOVASCULAR NEUROMODULATION

Global Headquarters

One St. Jude Medical Drive
St. Paul, Minnesota 55117
USA
+1 651 756 2000
+1 651 756 3301 Fax

Atrial Fibrillation Division

One St. Jude Medical Drive
St. Paul, Minnesota 55117
USA
+1 651 756 2000
+1 651 756 3301 Fax

SJM Coordination Center BVBA

The Corporate Village
Da Vincilaan 11 Box F1
1935 Zaventem
Belgium
+32 2 774 68 11
+32 2 772 83 84 Fax

St. Jude Medical Brasil Ltda.

Rua Itapeva, 538
5º ao 8º andar
01332-000 – São Paulo – SP
Brazil
+55 11 5080 5400
+55 11 5080 5423 Fax

St. Jude Medical (Hong Kong) Ltd

Suite 1608, 16/F Exchange Tower
33 Wang Chiu Road
Kowloon Bay, Kowloon
Hong Kong SAR
+852 2996 7688
+852 2956 0622 Fax

St. Jude Medical Japan Co., Ltd.

Shiodome City Center 15F
1-5-2, Higashi-Shinbashi
Minato-ku, Tokyo 105-7115
Japan
Tel: +81 3 6255 6370
Fax: +81 3 6255 6371

SJMprofessional.com



Rx Only

Brief Summary: Please review the Instructions for Use prior to using these devices for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

ViewMate, ViewFlex, Agilis, ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2012 St. Jude Medical. All Rights Reserved.

Item 100071833EN Rev A

